# ENFit<sup>®</sup> to Funnel Adaptor with cap

INSTRUCTIONS FOR USE

**GBUK Group Ltd.**  
Woodland House,  
Blackwood Hall Business Park,  
North Duffield, Selby,  
North Yorkshire, YO8 5DD, UK  
**EU: +44 (0)1757 282 945**  
**US: 18665 289 422**  
[www.gbukenteral.com](http://www.gbukenteral.com)  
[info@gbukenteral.com](mailto:info@gbukenteral.com)

**EC REP** **QualRep Services BV**  
Utrechtseweg 310 B42  
6812 AR Arnhem,  
The Netherlands

**CH REP** **Swiss AR Services AG**  
Industriestrasse 47  
6300 Zug, Switzerland

	<b>Page</b>
<b>EN</b> ISOSAF™ ENFit® to Funnel Adaptor with cap	4
<b>FR</b> Adaptateur ENFit® à entonnoir avec capuchon	5
<b>IT</b> Adattatore da ENFit® a imbuto con cappuccio	6
<b>ES</b> Adaptador ENFit® a embudo con tapa	7
<b>DE</b> ENFit® an Trichter-Adapter mit Verschlusskappe	8
<b>DA</b> ENFit® til tragt-adapter med hætte	9
<b>NO</b> ENFit® til traktadapteren med hette	10
<b>SV</b> ENFit®-till-tratt-adapter med lock	11
<b>NL</b> ENFit®-naar-trechteradapter met dop	12
<b>PL</b> Adapter ENFit® do lejka z zatyczką	13
<b>PT</b> Adaptador ENFit® para funil com tampa	14
<b>EL</b> ENFit® προς προσαρμογέα χοάνης με πώμα	15
<b>TR</b> ENFit®'ten Huniye kapaklı adaptör	16
<b>CS</b> Adaptér ENFit® k nálevce s víčkem	17
<b>ET</b> ENFit® korgiga lehteradapterile	18
<b>HU</b> ENFit®-tölcséres adapter, kupakkal	19
<b>LV</b> ENFit® piltuves adapteris ar vāciņu	20
<b>MT</b> ENFit® għal adapter ta' Lenbut bit-tapp	21
<b>RO</b> Adaptor ENFit® la pâlnie cu capac	22
<b>SL</b> Adapter ENFit® za lijak s pokrovčkom	23

EN

## ISOSAF™ ENFit® to Funnel Adaptor with cap

### Instruction For Use

Please read these instructions before use.

#### DEVICE DESCRIPTION

The UAD2ISO enteral transition adaptor is a single use medical device, provided non-sterile. The appropriate portion of the device is compliant with ISO 80369-3, allowing the connection between ISO 80369-3 compliant devices, and non-ISO 80369-3 compliant devices.

#### PATIENT POPULATION

Intended for patients requiring enteral administration of fluids, feed, or medication.

#### INTENDED USE

Intended to facilitate enteral specific connections between ISO 80369-3 compliant connectors and non-ISO 80369-3 compliant legacy enteral connectors.

#### INTENDED USER

The device is intended to be used by clinical personnel and others with appropriate training and knowledge.

#### INDICATION

The enteral transition adaptor is indicated for use with liquids and dispersible/soluble formulations.

#### CONTRAINDICATIONS

There are no known contraindications associated with the enteral transition adaptor.

#### WARNINGS & PRECAUTIONS

- This adaptor is intended only for connecting ENFit compatible accessories to non-ISO 80369-3 compliant legacy enteral connectors. The adaptor may only be used to attach ENFit compatible accessories to enteral feeding tubes.
- The transition adapter does not protect against inadvertent misconnections with connectors of other devices.
- Do not use with liquids or medicines that are not in dispersible/soluble formulations
- Routinely check for obstructions and blockages to ensure device performance. Discard the device if lack of performance, blockages or obstructions occur.
- Do not use if packaging damaged.
- The device is single use only. Do not reuse. Reuse may result in cross contamination and loss of device performance
- Do not reuse the enteral transition adaptor once fully disconnected from the enteral feeding system.
- If the adaptor is not properly pushed into the funnel connector, leakage may occur. During connection, apply suitable force to the adaptor to
- provide a secure connection. Do not push too hard, so that the funnel connector on the Enteral tube is damaged.
- If the adaptor is not properly locked via screw fit to an ISO 80369-3 Compliant female connector, leakage may occur. During connection, apply
- suitable force to the adaptor to provide a secure connection. Do not over tighten to prevent damage to the connectors.
- In the event of a serious incident involving the device the user should report the event to GBUK and the competent

authority in which the user/patient is established. In the event of a suspected device fault the user should report the event to GBUK, and the device should be returned to GBUK for examination.

#### PACK CONTENTS

1 x UAD2ISO Enteral Transition Adaptor

#### DIRECTIONS

- Inspect contents of packaging for damage.
- Connect the stepped conical portion of the adaptor into the Non-ISO 80369-3 compliant enteral connector on the Enteral tube via a compression/push fit.
- An ISO 80369-3 Compliant female connector may now be attached via a screw fit, to the Male ENFit connector of the UAD2ISO. Open the cap prior to use.
- Confirm any transitional connectors are connected to an enteral port and not an IV set. Once a secure connection is achieved, administration can commence.
- To close off the connection between ISO 80369-3 compliant connectors and non-ISO 80369-3 compliant legacy enteral connectors, unscrew the ISO 80369-3 Compliant male connector from the ISO 80369-3 Compliant female connector.
- Ensure the open port is closed off when the adaptor is not in use. Connect the tethered cap by screw fit to the ISO 80369-3 Compliant male connector.
- To resume administration, unscrew the tethered cap from the ISO 80369-3 Compliant male connector and refer to step 3.
- Clean the adaptor upon each assessment to ensure the external portion and the connector are free from debris. Cleaning should be as per that of the Enteral tube. Refer to the cleaning instructions provided with the Enteral tube.
- Inspect the UAD2ISO for visible signs of damage, leakage, and loss of device performance. If present, unscrew the ISO 80369-3 Compliant connection, and use moderate tension to detach the stepped conical portion of the adaptor from the non-ISO 80369-3 compliant enteral connector on the Enteral tube. Discard as per local hospital protocol.

#### CLEANING & DISINFECTION

Single use device. Dispose after use.

#### STORAGE CONDITIONS

Do not use if packaging is damaged. Keep out of direct sunlight. Keep dry.

#### DISPOSAL

Dispose of device in line with local policy.

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GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

GBUK Enteral is a trading name of **GBUK Enteral Ltd.**

**Rx Only** Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Symbol glossary available at [www.gbukenteral.com](http://www.gbukenteral.com).  
Not made with natural rubber latex.

ISOSAF™ is a trade mark of GBUK Group Ltd., registered in the UK and other countries and regions.  
ENFit® is a registered trade mark of GEDSA and used with their permission.

**GBUK Group Ltd.**

Woodland House,  
Blackwood Hall Business Park,  
North Duffield, Selby,  
North Yorkshire, YO8 5DD, UK  
**EU: +44 (0)1757 282 945**  
**US: 18665 289 422**  
[www.gbukenteral.com](http://www.gbukenteral.com)  
[info@gbukenteral.com](mailto:info@gbukenteral.com)

**QualRep Services BV**  
Utrechtseweg 310 B42  
6812 AR Arnhem,  
The Netherlands

**Swiss AR Services AG**  
Industriestrasse 47  
6300 Zug, Switzerland