GBUKBanana
Part of the GBUK Group ●●●●

Clip Slings

Single Patient Use

INSTRUCTIONS FOR USE

	Page
EN Single Patient Use Clip Slings	4
FR Harnais à clips à usage individuel	6
IT Imbracature con clip monouso	8
ES Eslinga con pinzas para uso en un solo paciente	10
DE Schlingen in Klammerform zur Verwendung bei einem einzelnen Patienten	12
DA Klemmeslynger – må kun anvendes til én patient.	14
NO Sele med klips for én pasient	16
SV Lyftsele med klämlås för enskild patient	18
NL Tilbanden met klemmen voor gebruik met één patiënt	20
PL Nosidła zatrzaskowe do użycia u jednego pacjenta	22
PT Cinto de transferência com clipe de fixação de utilização em paciente único	24
EL Αιώρα με κλιπ για χρήση σε μεμονωμένο ασθενή	26
TR Tek Hastada Kullanım İçin Klips Askı Bezleri	28
CS Závěs se sponou pro jednoho pacienta	30
ET Klambriga silmused ühe patsiendiga kasutamiseks	32
HU Egy beteg általi használatra való csíptetős hevederek	34
LV Individuālai lietošanai paredzēts slings ar fiksatoriem	36
MT Bragi għall-Klipps għal Użu fuq Pazjent Wieħed	38
RO Hamuri cu cleme pentru uz individual	40
SL Dvižna zanka s sponko za uporabo pri enem pacientu	42

GBUK Group Ltd.
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 288 587
US: 18665 289 422
www.gbukbanana.com
info@gbukbanana.com

EC REP **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP **Swiss AR Services AG**
Industriestrasse 47
6300 Zug, Switzerland

EN

Single Patient Use Clip Slings

Instructions for Use

Please read these instructions before use.

INTENDED USE

The GBUK Banana clip sling is intended to be used to transfer patients following a risk assessment. The sling should be used in conjunction with a mechanical lifting device (either mobile or ceiling hoist) to ensure safe lifting and repositioning of a patient.

INTENDED USER

The GBUK Banana clip sling is intended to be used by trained healthcare professionals or those with the appropriate training or knowledge. Refer to the Royal College of Nursing Guidelines.

INTENDED PATIENT POPULATION

Intended for patients who require lifting and/or repositioning following a suitable risk assessment.

INDICATIONS

The sling is indicated for use on patients with limited mobility requiring lifting and/or repositioning.

CONTRAINDICATIONS

If the sling is in contact with broken skin, do not use for longer than 30 days' to the lifetime of device. Do not use on patients whose weight exceeds the safe working load limit of the sling or hoist, whichever is lower.

WARNINGS & PRECAUTIONS

- Follow local manual handling policy from your organization, National Back Exchange or HSE guidelines.
- All handlers must be trained and competent when using the sling.
- Do not use if packaging is open or damaged.
- Do not use if the wash indicator has turned red, or if the slings show any signs of damage or wear.
- Do not launder or bleach the device.
- For single patient use only. Cross patient use of this device poses risk of infection and/or cross contamination.
- In the event of a serious incident involving the device the user should report the event to GBUK and the competent authority in which the user/patient is established. In the event of a suspected device fault the user should report the event to GBUK, and the device should be returned to GBUK for examination.

A risk assessment and visual inspection must be performed by competent person(s) before using the sling. This should take into account the following factors:

- The physical and medical status of the patient, including skin integrity,
- The status of the sling, including the appearance of the wash indicator and the presence of any damage, rips, tears, soil, loose stitching, or any other general signs of wear and tear,
- The compatibility of the sling, including if the correct size and type of the sling is being used,
- The capability of the caregivers using the equipment; and,
- The situation in which it is being used.

COMPATIBILITY

Please visit www.gbukgroup.com/clipslings for a list of compatible hoists.

SAFE WORKING LOAD LIMIT

The GBUK Banana clip slings have a safe working load of up to 230kg.

EXPECTED LIFETIME

LOLER regulations require you to carry out a visual inspection of the sling prior to each use. GBUK Banana recommends that you replace the sling after one month of use (assuming four transfers a day). However, the expected life of the sling is dependent on the actual conditions of use. If the sling shows any signs of damage or wear, dispose of it and select a new sling.

WASH INDICATOR

The sling may be wiped clean if required but should not be laundered. If the sling has been exposed to water, the wash indicator on the back of the sling will turn red and the sling must be disposed of.

MATERIALS

Polyester fabric, polyurethane fibre filling, nylon thread, nylon straps, nylon clips.

DIRECTIONS FOR USE

- Select the required sling from the box. Before removing the sling from its bag, place the sling against the middle of the patient's spine of the patient to verify that the correct size has been selected. The sling should cover the patient's crown to coccyx. Use the black stitching at the top and bottom of the sling as a visual guide. Width can also be measured at this point to check that it is suitable by placing the sling at the middle of the spine and then ensuring that this covers the width of the patient.
- Once the correct size has been selected, remove the sling from the bag and fit the sling by placing it behind the patient ensuring that the labels and positioning guides are facing away from the patient and are on the outside for the handler to see. Follow your local manual handling policy to position the patient within the sling.
- Ensure that the sling is fitted low enough down the patient's back so that the seam is just below the coccyx. The top seam of the sling should be around the crown of the head. A slide sheet may be used if required but should be removed once the sling is in the desired position.
- Once this is done, ensure the sling is centred at the back of the patient.
- Fit the leg sections by bringing each leg along each side and feed underneath the thighs. Ensure the material is flat and use a slide sheet if required.
- There are dignity straps on each leg. Feed a single leg strap through the dignity strap of the other leg if required. The patient handling risk assessment/care plan should reference whether this is required.
- Complete visual checks to make sure that all points above have been completed.
- Attach the shoulder and leg clips to the spreader bar of the hoist. Ensure that these clips are clicked securely into place.
- Before operating the hoist, ensure all clips have been secured to the spreader bars and that the patient's arms are inside the sling.
- On completion of the transfer, remove the sling from the patient following your local manual handling policy.

CLEANING & DISINFECTION

GBUK Banana clip slings are designed to be reused with the same patient i.e., the device is single patient use only. If lightly soiled, manually wipe down the sling as soon as possible using alcohol-based cleaning solutions (up to 70%) or normal disinfectant solutions until the sling is visibly clean. If the sling does not appear visibly clean and remains contaminated upon visual inspection, discard the sling in line with local hospital protocol.

STORAGE CONDITIONS

Storage of the sling should be in accordance with the local hospital policy and kept off the floor. Keep dry. Keep out of direct sunlight.

DISPOSAL

Dispose of device in line with local hospital policy.


Image 1

Not made with natural rubber latex.

Symbol glossary available at www.gbukbanana.com

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
Ni izdelano iz naravnega lateksa.

Slovar simbolov najdete na naslovu www.gbukbanana.com

GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

GBUK Banana is a trading name of **GBUK Group Ltd.**

Symbol glossary available at www.gbukbanana.com
Not made with natural rubber latex.

Caution: Federal law restricts this device to sale by or on the order of a manual handler, trained healthcare professional and others appropriately trained in use of the device.

GBUK Group Ltd.

Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 288 587
US: 18665 289 422
www.gbukbanana.com
info@gbukbanana.com

QualRep Services BV
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

Swiss AR Services AG
Industriestrasse 47
6300 Zug, Switzerland