**GBUK Banana**  
Part of the GBUK Group ●●●●

# Loop Slings

Single Patient Use

**INSTRUCTIONS FOR USE**

	<b>Page</b>
<b>EN</b> Single Patient Use Loop Slings	4
<b>FR</b> Harnais à boucles à usage individuel	6
<b>IT</b> Imbracature ad anello monouso	8
<b>ES</b> Eslinga continua con pinzas para uso en un solo paciente	10
<b>DE</b> Schlingen in Schlaufenform zur Verwendung bei einem einzelnen Patienten	12
<b>DA</b> Løkkesynger – må kun anvendes til én patient.	14
<b>NO</b> Selesløyfe for bruk av én pasient	16
<b>SV</b> Lyftselar med öglor för enskild patient	18
<b>NL</b> Tilbanden met lussen voor gebruik met één patiënt	20
<b>PL</b> Nosidła pętlowe do użycia u jednego pacjenta	22
<b>PT</b> Cinto com anéis de utilização em paciente único	24
<b>EL</b> Αιώρες με θηλιές για χρήση σε μεμονωμένο ασθενή	26
<b>TR</b> Tek Hastada Kullanım İçin Hasta Taşıma Halkalı Askı Bezleri	28
<b>CS</b> Závěs se smyčkou pro jednoho pacienta	30
<b>ET</b> Rõngassilmus ühe patsiendiga kasutamiseks	32
<b>HU</b> Egy beteg általi használatra való hurok hevederek	34
<b>LV</b> Individuālais lietošanai paredzēts cilpveida slings	36
<b>MT</b> Bragi għall-Klipps għal Użu fuq Pazjent Wiehed	38
<b>RO</b> Hamuri cu inel pentru uz individual	40
<b>SL</b> Dvižna zanka z ušesi za uporabo pri enem pacientu	42

**GBUK Group Ltd.**  
Woodland House,  
Blackwood Hall Business Park,  
North Duffield, Selby,  
North Yorkshire, YO8 5DD, UK  
**EU: +44 (0)1757 288 587**  
**US: 18665 289 422**  
[www.gbukbanana.com](http://www.gbukbanana.com)  
[info@gbukbanana.com](mailto:info@gbukbanana.com)

**EC REP** **QualRep Services BV**  
Utrechtseweg 310 B42  
6812 AR Arnhem,  
The Netherlands

**CH REP** **Swiss AR Services AG**  
Industriestrasse 47  
6300 Zug, Switzerland

EN

## Single Patient Use Loop Slings

### Instructions for Use

Please read these instructions before use.

#### INTENDED USE

The GBUK Banana loop sling is intended to be used to transfer patients following a risk assessment. The sling should be used in conjunction with a mechanical lifting device (either mobile or ceiling hoist) to ensure safe lifting and repositioning of a patient.

#### INTENDED USER

The GBUK Banana loop sling is intended to be used by trained healthcare professionals or those with the appropriate training or knowledge. Refer to the Royal College of Nursing Guidelines.

#### INTENDED PATIENT POPULATION

Intended for patients who require lifting and/or repositioning following a suitable risk assessment.

#### INDICATIONS

The sling is indicated for use on patients with limited mobility requiring lifting and/or repositioning.

#### CONTRAINDICATIONS

If the sling is in contact with broken skin, do not use for longer than 30 days' to the lifetime of device. Do not use on patients whose weight exceeds the safe working load limit of the sling or hoist, whichever is lower.

#### WARNINGS & PRECAUTIONS

- Follow local manual handling policy from your organization, National Back Exchange or HSE guidelines.
- All handlers must be trained and competent when using the sling.
- Do not use if packaging is open or damaged.
- Do not use if the wash indicator has turned red, or if the slings show any signs of damage or wear.
- Do not launder or bleach the device.
- For single patient use only. Cross patient use of this device poses risk of infection and/or cross contamination.
- In the event of a serious incident involving the device the user should report the event to GBUK and the competent authority in which the user/patient is established. In the event of a suspected device fault the user should report the event to GBUK, and the device should be returned to GBUK for examination

A risk assessment and visual inspection must be performed by competent person(s) before using the sling. This should take into account the following factors:

- The physical and medical status of the patient, including skin integrity,
- The status of the sling, including the appearance of the wash indicator and the presence of any damage, rips, tears, soil, loose stitching, or any other general signs of wear and tear,
- The compatibility of the sling, including if the correct size and type of the sling is being used,
- The capability of the caregivers using the equipment; and,
- The situation in which it is being used

#### COMPATABILITY

Slings that have loop attachments are suitable for use with a 2 or 4 point loop spreader bar.

#### SAFE WORKING LOAD LIMIT

The GBUK Banana loop slings have a safe working load of up to 400 kg for the XXL size, and 275 kg for all other sizes.

#### EXPECTED LIFETIME

LOLER regulations require you to carry out a visual inspection of the sling prior to each use. GBUK Banana recommends that you replace the sling after 6 months of use. However, the expected life of the sling is dependent on the actual conditions of use. If the sling shows any signs of damage or wear, dispose of it and select a new sling.

#### WASH INDICATOR

The sling may be wiped clean if required but should not be laundered. If the sling has been exposed to water, the wash indicator on the back of the sling will turn red and the sling must be disposed of.

#### MATERIALS

Polyester fabric, polyurethane fibre filling, nylon thread, nylon straps.

#### DIRECTIONS FOR USE

- Select the required sling from the box. Before removing the sling from its bag, place the sling against the middle of the patient's spine of the patient to verify that the correct size has been selected. The sling should cover the patient's crown to coccyx. Use the black stitching at the top and bottom of the sling as a visual guide. Width can also be measured at this point to check that it is suitable by placing the sling at the middle of the spine and then ensuring that this covers the width of the patient.
- Once the correct size has been selected, remove the sling from the bag and fit the sling by placing it behind the patient ensuring that the labels and positioning guides are facing away from the patient and are on the outside for the handler to see. Follow your local manual handling policy to position the patient within the sling.
- Ensure that the sling is fitted low enough down the patient's back so that the seam is just below the coccyx. The top seam of the sling should be around the crown of the head. A slide sheet may be used if required but should be removed once the sling is in the desired position.
- Once this is done, ensure the sling is centred at the back of the patient.
- Fit the leg sections by bringing each leg along each side and feed underneath the thighs. Ensure the material is flat and use a slide sheet if required.
- There are dignity straps on each leg. Feed a single leg strap through the dignity strap of the other leg if required. The patient handling risk assessment/care plan should reference whether this is required.
- Complete visual checks to make sure that all points above have been completed.
- Attach the shoulder and leg loops to the spreader bar of the hoist. The colour coded loops for each patient should be selected in accordance with the local manual handling policy and as per individual risk assessment/ care plan.
- Before operating the hoist, ensure all loops have been secured to the spreader bars and that the patient's arms are inside the sling.
- On completion of the transfer, remove the sling from the patient following your local manual handling policy.

#### CLEANING & DISINFECTION

GBUK Banana loop slings are designed to be reused with the same patient i.e., the device is single patient use only. If lightly soiled, manually wipe down the sling as soon as possible using alcohol-based cleaning solutions (up to 70%) or normal disinfectant solutions until the sling is visibly clean. If the sling does not appear visibly clean and remains contaminated upon visual inspection, discard the sling in line with local hospital protocol.

#### STORAGE CONDITIONS

Storage of the sling should be in accordance with the local hospital policy and kept off the floor. Keep dry. Keep out of direct sunlight.

#### DISPOSAL

Dispose of device in line with local hospital policy.

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
GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

GBUK Banana is a trading name of **GBUK Group Ltd.**

Symbol glossary available at [www.gbukbanana.com](http://www.gbukbanana.com)  
Not made with natural rubber latex.

**Caution: Federal law restricts this device to sale by or on the order of a manual handler, trained healthcare professional and others appropriately trained in use of the device.**

**GBUK Group Ltd.**

Woodland House,  
Blackwood Hall Business Park,  
North Duffield, Selby,  
North Yorkshire, YO8 5DD, UK  
**EU: +44 (0)1757 288 587**  
**US: 18665 289 422**  
[www.gbukbanana.com](http://www.gbukbanana.com)  
[info@gbukbanana.com](mailto:info@gbukbanana.com)

**QualRep Services BV**  
Utrechtseweg 310 B42  
6812 AR Arnhem,  
The Netherlands

**Swiss AR Services AG**  
Industriestrasse 47  
6300 Zug, Switzerland