LO-RAISER™

INSTRUCTIONS FOR USE

GBUK Group Ltd.
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 288 587
US: 18665 289 422
www.gbukbanana.com
info@gbukbanana.com

EC REP **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP **Swiss AR Services AG**
Industriestrasse 47,
6300 Zug,
Switzerland

	Page
EN LO-RAISER™	4
FR LO-RAISER™	6
IT LO-RAISER™	8
ES LO-RAISER™	10
DE LO-RAISER™	12
DA LO-RAISER™	14
NO LO-RAISER™	16
SV LO-RAISER™	18
NL LO-RAISER™	20
PL LO-RAISER™	22
PT LO-RAISER™	24
EL LO-RAISER™	26
TR LO-RAISER™	28
CS LO-RAISER™	30
ET LO-RAISER™	32
HU LO-RAISER™	34
LV LO-RAISER™	36
MT LO-RAISER™	38
RO LO-RAISER™	40
SL LO-RAISER™	42

EN

LO-RAISER™

Intended use

The Lo-Raiser is designed to support a person in moving from the floor to a seated position, or from a seated position to a standing position.

Indication for use

The Lo-Raiser is indicated for use on people who require assistance when moving from the floor to standing (100kg) or sitting to standing (200kg)

Description

Designed to help manually raise a person who has some lower body strength. The key to the use of the Lo-Raiser is to use weight transference to support the person up from the floor/chair/bed etc. The Lo-Raiser has a Safe Working Load (SWL) limit of 100kg from floor to standing and 200kg from sitting to standing.

Intended Target Group

This device is suitable for people who require assistance in raising from the floor to standing or sitting to standing. If the device is being used to raise a person from the floor, the person must be able to long sit with or without support and partially bend their knees. If the device is being used to raise a person from a seated position, the person must be able to take their own weight once standing.

Intended Users

The Lo-Raiser is intended to be used by trained healthcare professionals or those with the appropriate training or knowledge.

Contraindications

Do not use with a person who exceeds the SWL limit. The device is contraindicated for people with insufficient lower body mobility or strength. Do not use in direct contact with broken skin.

Warnings and Precautions

- The SWL limit is not a guide to the safe weight to handle. Assess the person's physical and medical status.
- Do not lift upwards with the straps or pull too rapidly.
- Follow local manual handling policy from your hospital, Royal College of Nursing or National Back Exchange guidelines.

A risk assessment and visual inspection must be performed by competent person(s) before using the Lo-Raiser. This should consider:

- The physical and medical status of the person, including skin integrity, weight, height, physical and cognitive ability.
- The status of the device. The material must be intact and not exhibit any splits or tears. The stitching on the handles should be inspected to ensure the stitching is intact. Do not use the device if it appears damaged, wet, or worn.
- The individual capability of the handlers and the environment. The footwear of the person and the surface should prevent slipping.

In the event of a serious incident involving the device, the user should report the event to GBUK and the competent authority in which the user/person is established. In the event of a suspected device fault, the user should report the event to GBUK, and the device should be returned to GBUK for examination.

DIRECTIONS FOR USE

FROM THE FLOOR

1. Slide the base of Lo-Raiser under the person's bottom. The person's bottom and base of the sacrum should be directly on the anti-slip area. You may need to rock or lean the person. You may also need a second handler to assist.
2. Fasten the belt strap and clip into place. Adjust to fit comfortably and securely. Pass the braces over the person's shoulders and cross the straps over the person's chest. The yellow clips should be connected to each other. The black clips should also be connected to each other. Tighten the straps if necessary while ensuring the person's comfort. Ensure that there is an equal length of handling strap/ladder on either side of the person.
3. The person will need to bend their knees slightly. Ensure the person's feet are flat on the floor before attempting to stand them. Dependent on the risk assessment, the feet can be blocked by the handler's foot.
4. Position yourself on either side of the person, adopting a 'lunge stance' with your outside leg forward, inside knee bent.
5. Your feet should be positioned around the area near the person's feet, and the outside foot should be near the person's hips/bottom. Adjust accordingly to accommodate your size.
6. Hold the ladder handles with your outside hand where you feel most comfortable. Apply moderate tension to the ladder before the stand process. If the person is tall, you may need to hold a handle closer to the hip. With the inside hand, hold the handle on the shoulder section of the braces. The person can hold onto your arms or any straps if they feel more comfortable. They can push down on the arms of a chair or on their own knees.
7. Co-ordinate the move using 'Ready, Steady, Stand/Step/Slide'. Move your outside leg back so that it is in line with your inside leg, levering the person from the floor/chair (Image 2). Your upper body becomes upright to help with leverage. The technique is a weight transference to lever the person, using the handler's weight with locked arms. Do not lift upwards with the straps or pull too rapidly. This will pull the Lo-Raiser out of position and you will have to start again.
8. The person can be assisted into a standing position, or, with the assistance of a 3rd or 4th handler, moved onto a chair/wheelchair. The 3rd handler pushes the chair underneath the person's bottom when they have been assisted partly up off the floor. The 4th handler may be holding the person's feet.

FROM A CHAIR

1. Apply the Lo-Raiser behind the person, ensuring that the lower edge of the Lo-Raiser is positioned slightly underneath the person's bottom, near the base of their sacrum. You may also need a second handler to assist.

2. Encourage the person to push down with their arms and instigate natural movement from sit to stand.
3. Use the same step motion to assist the person up, using 'ready, steady – stand' to guide movement (Image 3). Ensure the person's feet are flat on the floor before attempting to stand them. The Lo-Raiser can be used to 'walk' the person's bottom to the edge of the chair to aid in the movement.

Cleaning & disinfection

Suitable for laundering. If heavily soiled machine wash at 60 degrees Celsius as soon as possible. Leave to air dry. If lightly soiled manually wipe down the device until the device is visibly clean using alcohol-based cleaning solutions (up to 70%) or normal disinfectant solutions. If the device does not appear visibly clean, repeat the cleaning step. If the device remains contaminated upon visual inspection following multiple cleaning attempts, discard of the device in line with local hospital protocol. Do not bleach, dry clean or iron.

Storage and handling

Store the Lo-Raiser in the bag provided and keep off the floor. Keep dry and out of direct sunlight.

Disposal procedures

Follow your hospital protocol.

Not made with natural rubber latex.

Symbol glossary available at www.gbukbanana.com

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
slika 2

slika 3

GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

GBUK Banana is a trading name of **GBUK Group Ltd.**

Symbol glossary available at www.gbukbanana.com
Not made with natural rubber latex.

Caution: Federal law restricts this device to sale by or on the order of a manual handler, trained healthcare professional and others appropriately trained in use of the device.

GBUK Group Ltd.

Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 288 587
US: 18665 289 422
www.gbukbanana.com
info@gbukbanana.com

QualRep Services BV
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

Swiss AR Services AG
Industriestrasse 47
6300 Zug, Switzerland