Flo™ Air Transfer Device

Single Patient Use

INSTRUCTIONS FOR USE

	Page
EN GBUK Banana FLO™ Air Transfer Device	4
FR Dispositif de transfert pneumatique GBUK Banana FLO™	6
IT Dispositivo ad aria Banana FLO™ GBUK per il trasferimento dei pazienti	8
ES Producto sanitario neumático para tareas de transferencia Banana FLO™ GBUK	10
DE Luftunterstützte GBUK Banana FLO™-Vorrichtung für die Patientenverlegung	12
DA GBUK Banana FLO™ Air Transfer Device	14
NO Luftoverføringsenhet fra GBUK Banana FLO™	16
SV GBUK Banana FLO™ Luftöverföringsenhet	18
NL GBUK Banana FLO™-pneumatisch transferhulpmiddel	20
PL Urządzenie pneumatyczne do transferu GBUK Banana FLO™	22
PT Dispositivo de transferência pneumático GBUK Banana FLO™	24
EL Συσκευή μεταφοράς αέρα GBUK Banana FLO™	26
TR GBUK Banana FLO™ Şişme Transfer Cihazı	28
CS Prostředek pro přenos vzduchem GBUK Banana FLO™	30
ET GBUK Banana FLO™ täispuhutav teilsaldusvahend	32
HU GBUK Banana FLO™ Beteg mozgató csúszópád	34
LV GBUK Banana FLO™ pneimatiska pārvietošanas ierīce	36
MT Tagħmir għat-Trasferiment tal-Arja FLO™ tal-GBUK Banana	38
RO Dispozitiv de transfer de aer GBUK Banana FLO™	40
SL GBUK Banana FLO™ Air naprava za prenašanje	42

GBUK Group Ltd.
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 288 587
US: 18665 289 422
www.gbukbanana.com
info@gbukbanana.com

QualRep Services BV
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

Swiss AR Services AG
Industriestrasse 47
6300 Zug, Switzerland

EN

GBUK Banana FLO™ Air Transfer Device

Intended use

These devices are intended to assist in laterally transferring patients from one surface to another.

Indications for use

This equipment is indicated for patients within the patient target group.

- Patients who are dependent, unable to, or not required to participate in their transfer.
- Patients whose body weight and size pose a significant risk or care management issue to the patient or staff during performance of routine nursing care.

Description

These devices have been designed to reduce risks and forces and make patient transfers easier and safer. These devices have a safe working load (SWL) limit of 453 kg (1,000 lbs).

Product Code	Product Description
B-FLO34	Banana FLO™ 34" Single Patient Use Air Transfer Device
B-FLO39	Banana FLO™ 39" Single Patient Use Air Transfer Device
B-FLO50	Banana FLO™ 50" Single Patient Use Air Transfer Device
B-FLOSHORT34	Banana FLO™ 34" Short Single Patient Use Air Transfer Device

Patient target group

These devices are suitable for patients who meet the indications for use.

Intended users

These devices are intended to be used by trained healthcare professionals or those with the appropriate training or knowledge.

Contra-indications

These devices are contraindicated for:

- Patients with thoracic, cervical or lumbar fractures
- Patients whose body shapes do not fit within the physical boundaries of the mattress

Do not use in direct contact with broken skin.

Do not use with patients who exceed the SWL limit.

Warnings and precautions

- Follow local manual handling policy from your hospital or Royal College of Nursing guidelines.
- Do not use the device for lifting.
- Do not exceed the 453 kg /1000 lb SWL limit.
- For single patient use only. Cross patient use of this device poses a risk of infection and/or cross contamination.
- Do not use device to bridge a gap larger than 7.6 cm / 3".
- Ensure that all caster brakes (if any) are engaged prior to transfer.

- Do not leave a patient unattended on an inflated transfer device.
- Do not attempt to move a patient on an un-inflated transfer device.
- Ensure that the side-rail on the opposite side of the receiving surface is in the upright and locked position before transfer.
- Ensure that all patient support systems such as I.V. lines or oxygen hoses are free to transfer with the patient.
- Always use at least two people during patient transfers.
- Ensure that the patient is placed on the non-perforated side of the air transfer device with feet at the label end of the mattress.

A risk assessment and visual inspection must be performed by competent person(s) before using the Air Transfer Device. This should consider:

- The physical and medical status of the patient, including skin integrity.
- The status of the device. The material must be dry, intact, and not exhibit any splits or tears. The stitching on the handles should be inspected to ensure the stitching is intact. Do not use the device if it appears damaged, wet, or worn.
- The suitability of the device, including size.

In the event of a serious incident involving the device, the user should report the event to GBUK and the competent authority in which the user/patient is established. In the event of a suspected device fault, the user should report the event to GBUK, and the device should be returned to GBUK for examination.

Operating instructions

The device is compatible with air pumps with an air supply of $\geq 2,500$ L/min.

1. Follow your local manual handling policy for manoeuvres listed below.
2. Ensure that any brakes are engaged. If possible, the side rail on the opposite side of the receiving surface should be raised and locked in position. It is advantageous to have the surface you are transferring to lower than the surface you are transferring from.
3. Position the device beneath the patient.
4. Gently secure patient straps across the patient.
5. Check there is nothing in the transfer area that will obstruct or damage the device.
6. Check all patient support systems such as I.V. lines or oxygen hoses are free to transfer with the patient.
7. Connect the air supply hose to either inlet near the patient's feet, ensuring that the hose is not kinked and remains unobstructed. Snap the hose into place.
8. Turn on air supply and allow device to inflate fully. When using a 50" device, two air supplies may be connected simultaneously.
9. Transfer the patient. Diagonal transfers are usually preferred as it reduces the effort required and provides a sense of security to the patient.
10. Disconnect the hose and remove the device from underneath the patient.

Storage and handling

Keep dry and out of direct sunlight.

Disposal procedures

Follow your hospital protocol.

Cleaning

These devices should not be laundered. If lightly soiled, manually wipe down the device as soon as possible using alcohol-based cleaning solutions (up to 70%) or normal disinfectant solutions until the device is visibly clean. If the device does not appear visibly clean and remains contaminated upon visual inspection, discard the device in line with local hospital protocol.

Not made with natural rubber latex.

Symbol glossary available at www.gbukbanana.com

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Ni izdelano iz naravnega lateksa.

Slovar simbolov najdete na naslovu www.gbukbanana.com

GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

GBUK Banana is a trading name of **GBUK Group Ltd.**

Symbol glossary available at www.gbukbanana.com
Not made with natural rubber latex.

Flo is a trade mark of GBUK Group Ltd., registered in the UK and other countries and regions.

GBUK Group Ltd.

Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0)1757 288 587
US: 18665 289 422
www.gbukbanana.com
info@gbukbanana.com

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