**GBUK Healthcare**  
Part of the GBUK Group ●●●●

DASH 6™

Drawing Up

Filter Straw (BNS)

**INSTRUCTIONS FOR USE**

 **GBUK Group Ltd.**  
Woodland House,  
Blackwood Hall Business Park,  
North Duffield, Selby,  
North Yorkshire, YO8 5DD, UK  
**EU: +44 (0) 1757 288 587**  
**US: 18665 289 422**  
[www.gbukhealthcare.com](http://www.gbukhealthcare.com)  
[info@gbukhealthcare.com](mailto:info@gbukhealthcare.com)

**EC REP** **QualRep Services BV**  
Utrechtseweg 310 B42  
6812 AR Arnhem,  
The Netherlands

**CH REP** **Swiss AR Services GmbH**  
Industriestrasse 47,  
6300 Zug,  
Switzerland

	<b>Page</b>
<b>EN</b> DASH 6™ Drawing Up Filter Straw	4
<b>FR</b> Paille De Prélèvement Avec Filtre DASH 6™	5
<b>IT</b> Cannula Filtrante Per Aspirazione DASH 6™	6
<b>ES</b> Cánula De Extracción Con Filtro DASH 6™	7
<b>DE</b> DASH 6™ Aufziehfilterhalm	8
<b>DA</b> DASH 6™ Slange Til Optrækning Med Filter	9
<b>NO</b> DASH 6™ Oppsugingsrør Med Filter	10
<b>SE</b> DASH 6™ Uppsugningrör Med Filter	11
<b>NL</b> DASH 6™-Aanzuigfilterslang	12
<b>PL</b> Rurka Do Pobierania Z Filtrem DASH 6™	13
<b>PT</b> Canudo Do Filtro De Extração DASH 6™	14
<b>EL</b> Σωληναριο Φιλτρου Αντλησης DASH 6™	15
<b>TR</b> DASH 6™ Çekme Filtresi Pip	16
<b>CS</b> DASH 6™ – Filtrační Rourka Pro Nasávání	17
<b>ET</b> DASH 6™ Sissetõmbe Kõrsfilter	18
<b>HU</b> DASH 6™ Szűrős Felszívó Szál	19
<b>LV</b> DASH 6™ Iesūkšanas Salmiņš Ar Filtru	20
<b>MT</b> Strò Ta' Filtru Li Jigbed 'Il Fuq DASH 6™	21
<b>RO</b> DASH 6™ Pai Filtru De Extragere	22
<b>SL</b> Filtrska Cevka DASH 6™ Za Zajemanje	23
<b>SK</b> Filtračná Trubica DASH 6™	24
<b>BG</b> Филтърна Сламка За Изтегляне DASH 6™	25
<b>LT</b> DASH 6™ Pritraukimo Šiaudelis Su Filtru	26

EN

## DASH 6™ DRAWING UP FILTER STRAW

### Instructions For Use

Please read these instructions before use.

### INTENDED USE

The DASH 6™ Drawing Up Filter Straw is intended to be used with ISO 80369-6 NRFit® neuraxial syringes for the drawing up of neuraxial medication and anaesthetic.

### INDICATION

The device is used with an ISO 80369-6 NRFit® compliant neuraxial syringe and is used for drawing up and filtering neuraxial medication and anaesthetic for subsequent administration to a patient.

### INTENDED USER

The device is intended to be used by a trained healthcare professional.

### CONTRAINDICATIONS

The device is not intended for parenteral or enteral use.

### WARNINGS & PRECAUTIONS

- Do not use if packaging is damaged.
- Do not re-use or reprocess the device as this may lead to leakage or infection.
- Do not use if the tip is visibly blocked or if the device is damaged.
- The device is single use only. Reuse may result in loss of device performance and/or infection.

### DIRECTIONS


Always follow the hospital protocol when using the device. The device is to be used in conjunction with the ISO 80369-6 NRFit® syringe. Please follow the instructions for use as per the ISO 80369-6 NRFit® syringe provided for drawing up of neuraxial medication and anaesthetic. Attach the device to an ISO 80369-6 NRFit® syringe by twisting the hub of the Drawing Up Filter Straw onto the ISO 80369-6 NRFit® compliant syringe connector.

### STORAGE CONDITION

Keep dry and away from direct sunlight.

### DISPOSAL

Dispose of the device as per hospital protocol.

### AVAILABLE PRODUCTS

#### DESCRIPTION:

- DASH 6™ ISO 80369-6 NRFit® Compliant Drawing Up Filter Straw 5µ x 100mm

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GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

\*Trademark or Registered Trademark of GBUK Group Ltd.  
[www.gbukgroup.com](http://www.gbukgroup.com)

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Symbol glossary available at [www.gbukhealthcare.com](http://www.gbukhealthcare.com).

**GBUK Group Ltd.**

Woodland House,  
Blackwood Hall Business Park,  
North Duffield, Selby,  
North Yorkshire, YO8 5DD, UK  
**EU: +44 (0) 1757 288 587**  
**US: 18665 289 422**  
[www.gbukhealthcare.com](http://www.gbukhealthcare.com)  
[info@gbukhealthcare.com](mailto:info@gbukhealthcare.com)

**QualRep Services BV**  
Utrechtseweg 310 B42  
6812 AR Arnhem,  
The Netherlands

**Swiss AR Services GmbH**  
Industriestrasse 47,  
6300 Zug,  
Switzerland