DASH 6™
Epidural Catheter
Feeder

INSTRUCTIONS FOR USE

GBUK Group Ltd.
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0) 1757 288 587
US: 18665 289 422
www.gbukhealthcare.com
info@gbukhealthcare.com

EC REP **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP **Swiss AR Services AG**
Industriestrasse 47,
6300 Zug,
Switzerland

	Page
EN DASH 6™ Epidural Catheter Feeder	4
FR Distributeur De Cathéter Épidural DASH 6™	5
IT Alimentatore Per Catetere Epidurale DASH 6™	6
ES Guía De Catéter Epidural DASH 6™	7
DE DASH 6™ Epiduralkatheterzuführung	8
DA DASH 6™ Epidural Kateterindføringsanordning	9
NO DASH 6™ Epiduralkatetermater	10
SE DASH 6™ Epiduralkatetermatrare	11
NL DASH 6™-Epidurale Kathetergeleider	12
PL Prowadnica Do Cewnika Zewnařrzoponowego DASH 6™	13
PT Alimentador De Cateter Epidural DASH 6™	14
EL Τροφοδοτης Επισκληριδιου Καθετηρα DASH 6™	15
TR DASH 6™ Epidural Kateter Besleyi	16
CS DASH 6™ – Zavaděč Epidurálního Katétru	17
ET DASH 6™ Epiduraalkateetri Annusti	18
HU DASH 6™ Epidurális Katéterbevezető	19
LV DASH 6™ Epidurālā Katetra Padevējs	20
MT Feeder Tal-Kateter Epidurali DASH 6™	21
RO Alimentator Cu Cateter Epidural DASH 6™	22
SL Napeljevalnik Epiduralnega Katetra DASH 6™	23
SK Zavadzač Epidurálneho Katétra DASH 6™	24
BG Подаващо Устройство За Епидурален Катетър DASH 6™	25
LT DASH 6™ Epidurinio Katetrio TiekTuvas	26

EN

DASH 6™ EPIDURAL CATHETER FEEDER

Instructions For Use

Please read these instructions before use.

INTENDED USE

The DASH 6™ Epidural Catheter Feeder is intended to be used with ISO 80369-6 NRFit® neuraxial compliant devices to assist the insertion of an epidural needle into the epidural space.

INDICATION

The component is used with an ISO 80369-6 NRFit® compliant epidural needle and epidural catheter and is used to assist the catheter placement into the epidural space.

INTENDED USER

The device is intended to be used by a trained healthcare professional.

CONTRAINDICATIONS

The device is not intended for parenteral or enteral use.

WARNINGS & PRECAUTIONS

- Do not use if packaging is damaged.
- The component is provided non-sterile and it is the responsibility of the kit manufacturer to ensure the component can be packaged and sterilised under their internal processes.
- It is the responsibility of the kit manufacturer to ensure that the component is compatible with the devices in their kits.
- The device is for single use only. Reuse may result in loss of device performance and/or infection.

DIRECTIONS

Always follow the hospital protocol when using the device. The device is to be used in conjunction with the ISO 80369-6 NRFit® Epidural Needle and Catheter. Please follow the instructions for use as per the hospital procedure/instructions when applying the epidural needle/epidural catheter into the patient.


1. The Epidural Catheter is fed through the top of the feeder until the tip of the catheter can be seen on the ISO 80369-6 NRFit® connector tip.
2. The Epidural Catheter Feeder is then attached to the ISO 80369-6 NRFit® compliant Epidural Needle and the catheter is fed through the epidural needle to the necessary length.

STORAGE CONDITION

Keep dry and away from direct sunlight.

DISPOSAL

Dispose of the device as per hospital protocol.

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GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

[™]Trademark or Registered Trademark of GBUK Group Ltd.

www.gbukgroup.com

Rx Only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Symbol glossary available at www.gbukhealthcare.com.

GBUK Group Ltd.

Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0) 1757 288 587
US: 18665 289 422
www.gbukhealthcare.com
info@gbukhealthcare.com

QualRep Services BV
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

Swiss AR Services AG
Industriestrasse 47,
6300 Zug,
Switzerland