GBUKHealthcare
Part of the GBUK Group ●●●●

DASH 6™

Blunt Drawing Up
Filter Needle (BNS)

INSTRUCTIONS FOR USE

GBUK Group Ltd.
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0) 1757 288 587
US: 18665 289 422
www.gbukhealthcare.com
info@gbukhealthcare.com

EC REP **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP **Swiss AR Services GmbH**
Industriestrasse 47,
6300 Zug,
Switzerland

	Page
EN DASH 6™ Blunt Drawing Up Filter Needle	4
FR Aiguille Émoussée De Prélèvement Avec Filtre DASH 6™	5
IT Ago Filtrante Smussato Per Aspirazione DASH 6™	6
ES Aguja Roma De Extracción Con Filtro DASH 6™	7
DE DASH 6™ Stumpfe Ansaugfilterkanüle	8
DA DASH 6™ Stump Nål Til Optrækning Med Filter	9
NO DASH 6™ Butt Oppsugingskanyle Med Filter	10
SE Topa Filtriska Igla DASH 6™ Za Zajemanje	11
NL DASH 6™ Budget Suction Filter Needle	12
PL Agulha Do Filtro De Extração Com Ponta Romba DASH 6™	13
PT Agulha Do Filtro De Extração Com Ponta Romba DASH 6™	14
EL Βελονα Αντλησης Φιλτρου Με Αμβλεια Ακρη DASH 6™	15
TR DASH 6™ Künt Çekme Filtresi İğne	16
CS DASH 6™ – Tupá Filtrační Jehla Pro Nasávání	17
ET DASH 6™ Nüri Täitmise Filternõel	18
HU DASH 6™ Tompa, Szűrős Felszívótű	19
LV DASH 6™ Neasā Ievilkšanas Adata Ar Filtru	20
MT Labra Spuntata Ta' Filtru Li Jigbed 'Il Fuq DASH 6™	21
RO DASH 6™ Ac Bont Filtru De Extragere	22
SL Topa Filtriska Igla DASH 6™ Za Zajemanje	23
SK Tupá Nasávacia Ihla S Filtrom DASH 6™	24
BG Филтърна Игла Със Затъпен Врѝх За Изтегляне DASH 6™	25
LT DASH 6™ Buka Pritraukimo Adata Su Filtru	26

EN

DASH 6™ BLUNT DRAWING UP FILTER NEEDLE

Instructions For Use

Please read these instructions before use.

INTENDED USE

The DASH 6™ Blunt Drawing Up Filter Needle is intended to be used with ISO 80369-6 NRFit® neuraxial syringes for the drawing up of neuraxial medication and anaesthetic.

INDICATION

The device is used with an ISO 80369-6 NRFit® compliant neuraxial syringe and is used for drawing up and filtering neuraxial medication and anaesthetic for subsequent administration to a patient.

INTENDED USER

The device is intended to be used by a trained healthcare professional.

CONTRAINDICATIONS

The device is not intended for parenteral or enteral use.

WARNINGS & PRECAUTIONS

Do not use if packaging is damaged.

Do not reuse or reprocess the device.

Do not use if the tip is visibly blocked or if the device is damaged.

The device is single use only. Reuse may result in loss of device performance and/or infection.

DIRECTIONS

Always follow the hospital protocol when using the device.

The device is to be used in conjunction with the ISO 80369-6 NRFit® syringe. Please follow the instructions for use as per the ISO 80369-6 NRFit® syringe provided for drawing up of neuraxial medication and anaesthetic. Attach the device to an ISO 80369-6 NRFit® syringe by twisting the hub of the blunt drawing up filter needle onto the ISO 80369-6 NRFit® compliant syringe connector.

STORAGE CONDITION

Keep dry and away from direct sunlight.

DISPOSAL

Dispose of the device as per hospital protocol.

AVAILABLE PRODUCTS

DESCRIPTION:

- DASH 6™ ISO 80369-6 NRFit® Compliant Blunt Drawing Up Filter Needle 5µm, 18G x 50mm

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GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

*Trademark or Registered Trademark of GBUK Group Ltd.
www.gbukgroup.com

Rx Only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Symbol glossary available at www.gbukhealthcare.com.

 **GBUK Group Ltd.**
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0) 1757 288 587
US: 18665 289 422
www.gbukhealthcare.com
info@gbukhealthcare.com

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