**GBUK Healthcare**  
Part of the GBUK Group ●●●●

DASH 6™

Tuohy Borst

Adaptor

**INSTRUCTIONS FOR USE**

**GBUK Group Ltd.**  
Woodland House,  
Blackwood Hall Business Park,  
North Duffield, Selby,  
North Yorkshire, YO8 5DD, UK  
**EU: +44 (0) 1757 288 587**  
**US: 18665 289 422**  
[www.gbukhealthcare.com](http://www.gbukhealthcare.com)  
[info@gbukhealthcare.com](mailto:info@gbukhealthcare.com)

**EC REP** **QualRep Services BV**  
Utrechtseweg 310 B42  
6812 AR Arnhem,  
The Netherlands

**CH REP** **Swiss AR Services GmbH**  
Industriestrasse 47,  
6300 Zug,  
Switzerland

	<b>Page</b>
<b>EN</b> DASH 6™ Tuohy Borst Adaptor	4
<b>FR</b> Adaptateur DASH 6™ Tuohy Borst	5
<b>IT</b> Adattatore Tuohy Borst DASH 6™	6
<b>ES</b> Adaptador Tuohy Borst DASH 6™	7
<b>DE</b> DASH 6™ Tuohy Borst-Adapter	8
<b>DA</b> DASH 6™ Tuohy Borst Adapter	9
<b>NO</b> DASH 6™ Tuohy Borst-Adapter	10
<b>SE</b> Adapter Tuohy Borst DASH 6™	11
<b>NL</b> DASH 6™-Tuohy-Borstadapter	12
<b>PL</b> Adapter Tuohy Borst DASH 6™	13
<b>PT</b> Adaptador DASH 6™ Tuohy Borst	14
<b>EL</b> Προσαρμογέας DASH 6™ Tuohy Borst	15
<b>TR</b> DASH 6™ Tuohy Borst Adaptörü	16
<b>CS</b> DASH 6™ – Adaptér Tuohy-Borst	17
<b>ET</b> DASH 6™ Tuohy Borsti Adapter	18
<b>HU</b> DASH 6™ Tuohy-Borst Adapter	19
<b>LV</b> DASH 6™ Tuohy Borst Adapteris	20
<b>MT</b> Adapter Tuohy Borst DASH 6™	21
<b>RO</b> DASH 6™ Adaptor Tuohy Borst	22
<b>SL</b> Adapter Tuohy Borst DASH 6™	23
<b>SK</b> Adaptér DASH 6™ Tuohy Borst	24
<b>BG</b> Адаптер Tuohy Borst DASH 6™	25
<b>LT</b> DASH 6™ Tuohy Borst Adapteris	26

EN

## DASH 6™ TUOHY BORST ADAPTOR

### Instructions For Use

Please read these instructions before use.

### INTENDED USE

The DASH 6™ Tuohy Borst Adaptor is intended to be used with an epidural catheter to provide an ISO 80369-6 NRFit® compliant connection.

### INDICATION

The component is used with an ISO 80369-6 NRFit® compliant epidural catheter and is used to form an NRFit® connection, which will allow administration/infusion of neuraxial medication.

### INTENDED USER

The device is intended to be used by a trained healthcare professional.

### CONTRAINDICATIONS

The device is not intended for parenteral or enteral use.

### WARNINGS & PRECAUTIONS

- Do not use if packaging is damaged.
- The component is provided non-sterile and it is the responsibility of the kit manufacturer to ensure the component can be packaged and sterilised under their internal processes.
- It is the responsibility of the kit manufacturer to ensure that the component is compatible with the devices in their kits.
- The device is for single use only. Reuse may result in loss of device performance and/or infection.

### DIRECTIONS

Always follow the hospital protocol when using the device. The device is to be used in conjunction with the ISO 80369-6 NRFit® catheter. Please follow the instructions for use as per the hospital procedure/instructions when applying the epidural catheter into the patient.


1. Ensure the Tuohy Borst Adaptor is unscrewed. Feed the catheter through the small hole in Tuohy Borst Adaptor until the catheter cannot be fed in any further.
2. Secure the catheter by twisting the Tuohy Borst Adaptor until it is locked.

### STORAGE CONDITION

Keep dry and away from direct sunlight.

### DISPOSAL

Dispose of the device as per hospital protocol.

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GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

\*Trademark or Registered Trademark of GBUK Group Ltd.  
[www.gbukgroup.com](http://www.gbukgroup.com)

**Rx Only** Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Symbol glossary available at [www.gbukhealthcare.com](http://www.gbukhealthcare.com).

 **GBUK Group Ltd.**  
Woodland House,  
Blackwood Hall Business Park,  
North Duffield, Selby,  
North Yorkshire, YO8 5DD, UK  
**EU: +44 (0) 1757 288 587**  
**US: 18665 289 422**  
[www.gbukhealthcare.com](http://www.gbukhealthcare.com)  
[info@gbukhealthcare.com](mailto:info@gbukhealthcare.com)

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