GBUK Healthcare
Part of the GBUK Group ●●●●

DASH 6™

Tuohy Borst

Adaptor

INSTRUCTIONS FOR USE

GBUK Group Ltd.
 Woodland House,
 Blackwood Hall Business Park,
 North Duffield, Selby,
 North Yorkshire, YO8 5DD, UK
EU: +44 (0) 1757 288 587
US: 18665 289 422
www.gbukhealthcare.com
info@gbukhealthcare.com

EC REP **QualRep Services BV**
 Utrechtseweg 310 B42
 6812 AR Arnhem,
 The Netherlands

CH REP **Swiss AR Services AG**
 Industriestrasse 47,
 6300 Zug,
 Switzerland

	Page
EN DASH 6™ Tuohy Borst Adaptor	4
FR Adaptateur DASH 6™ Tuohy Borst	5
IT Adattatore Tuohy Borst DASH 6™	6
ES Adaptador Tuohy Borst DASH 6™	7
DE DASH 6™ Tuohy Borst-Adapter	8
DA DASH 6™ Tuohy Borst Adapter	9
NO DASH 6™ Tuohy Borst-Adapter	10
SE Adapter Tuohy Borst DASH 6™	11
NL DASH 6™-Tuohy-Borstadapter	12
PL Adapter Tuohy Borst DASH 6™	13
PT Adaptador DASH 6™ Tuohy Borst	14
EL Προσαρμογέας DASH 6™ Tuohy Borst	15
TR DASH 6™ Tuohy Borst Adaptörü	16
CS DASH 6™ – Adaptér Tuohy-Borst	17
ET DASH 6™ Tuohy Borsti Adapter	18
HU DASH 6™ Tuohy-Borst Adapter	19
LV DASH 6™ Tuohy Borst Adapteris	20
MT Adapter Tuohy Borst DASH 6™	21
RO DASH 6™ Adaptor Tuohy Borst	22
SL Adapter Tuohy Borst DASH 6™	23
SK Adaptér DASH 6™ Tuohy Borst	24
BG Адаптер Tuohy Borst DASH 6™	25
LT DASH 6™ Tuohy Borst Adapteris	26

EN

DASH 6™ TUOHY BORST ADAPTOR

Instructions For Use

Please read these instructions before use.

INTENDED USE

The DASH 6™ Tuohy Borst Adaptor is intended to be used with an epidural catheter to provide an ISO 80369-6 NRFit® compliant connection.

INDICATION

The component is used with an ISO 80369-6 NRFit® compliant epidural catheter and is used to form an NRFit® connection, which will allow administration/infusion of neuraxial medication.

INTENDED USER

The device is intended to be used by a trained healthcare professional.

CONTRAINDICATIONS

The device is not intended for parenteral or enteral use.

WARNINGS & PRECAUTIONS

- Do not use if packaging is damaged.
- The component is provided non-sterile and it is the responsibility of the kit manufacturer to ensure the component can be packaged and sterilised under their internal processes.
- It is the responsibility of the kit manufacturer to ensure that the component is compatible with the devices in their kits.
- The device is for single use only. Reuse may result in loss of device performance and/or infection.

DIRECTIONS

Always follow the hospital protocol when using the device. The device is to be used in conjunction with the ISO 80369-6 NRFit® catheter. Please follow the instructions for use as per the hospital procedure/instructions when applying the epidural catheter into the patient.


1. Ensure the Tuohy Borst Adaptor is unscrewed. Feed the catheter through the small hole in Tuohy Borst Adaptor until the catheter cannot be fed in any further.
2. Secure the catheter by twisting the Tuohy Borst Adaptor until it is locked.

STORAGE CONDITION

Keep dry and away from direct sunlight.

DISPOSAL

Dispose of the device as per hospital protocol.

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Notes

GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

[™]Trademark or Registered Trademark of GBUK Group Ltd.

www.gbukgroup.com

Rx Only Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Symbol glossary available at www.gbukhealthcare.com.

 **GBUK Group Ltd.**
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0) 1757 288 587
US: 18665 289 422
www.gbukhealthcare.com
info@gbukhealthcare.com

EC REP **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP **Swiss AR Services AG**
Industriestrasse 47,
6300 Zug,
Switzerland