GBUK Healthcare
Part of the GBUK Group ●●●●

DASH 6™

Syringe Cap

(BNS)

INSTRUCTIONS FOR USE

 **GBUK Group Ltd.**
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0) 1757 288 587
US: 18665 289 422
www.gbukhealthcare.com
info@gbukhealthcare.com

EC REP **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP **Swiss AR Services AG**
Industriestrasse 47,
6300 Zug,
Switzerland

	Page
EN DASH 6™ Syringe Cap	4
FR Capuchon De Seringue DASH 6™	5
IT Cappuccio Per Siringa DASH 6™	6
ES Tapa De Jeringa DASH 6™	7
DE DASH 6™ Spritzenverschlusskappe	8
DA DASH 6™ Sprøjtehætte	9
NO DASH 6™ Sprøytehette	10
SV DASH 6™ Sprutlock	11
NL DASH 6™-Afsluitdop	12
PL Zatyczka Do Strzykawki DASH 6™	13
PT Tampa De Seringa DASH 6™	14
EL Καπάκι Συριγγας DASH 6™	15
TR DASH 6™ Şiringa Kapağı	16
CS DASH 6™ Krytka Stríkačky	17
ET DASH 6™ Süstlakork	18
HU DASH 6™ Fecskendőkupak	19
LV DASH 6™ Šļircis Vāciņš	20
MT Tapp Tas-Siringa DASH 6™	21
RO DASH 6™ Capac Seringă	22
SL Pokrovček Brizgalke DASH 6™	23
SK Krytka Injekčnej Striekačky DASH 6™	24
BG Капачка За Спринцовка DASH 6™	25
LT DASH 6™ Švirkšto Gaubtelis	26

EN

DASH 6™ SYRINGE CAP

Instructions For Use

Please read these instructions before use.

INTENDED USE

The DASH 6™ Syringe Cap is intended to be used with ISO 80369-6 NRFit® neuraxial syringes for sealing the tip of the syringe.

INDICATION

The device is used with an ISO 80369-6 NRFit® compliant neuraxial syringe and is used to seal the contents of a syringe which is not intended for immediate use.

INTENDED USER

The device is intended to be used by a trained healthcare professional.

CONTRAINDICATIONS

The device is not intended for parenteral or enteral use.

WARNINGS & PRECAUTIONS

- Do not use if packaging is damaged.
- Do not re-use or reprocess the device as this may lead to leakage.
- The device is single use only. Reuse may result in loss of device performance and/or infection.

DIRECTIONS

Always follow the hospital protocol when using the device.

The device is to be used in conjunction with the ISO 80369-6 NRFit® syringe. Please follow the instructions for use as per the ISO 80369-6 NRFit® syringe provided for sealing the tip of the syringe.


Attach the device to an ISO 80369-6 NRFit® syringe by twisting the cap onto the ISO 80369-6 NRFit® compliant syringe connector.

STORAGE CONDITION

Keep dry and away from direct sunlight.

DISPOSAL

Dispose of the device as per hospital protocol.

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
Notes

GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

[™]Trademark or Registered Trademark of GBUK Group Ltd.

www.gbukgroup.com

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Symbol glossary available at www.gbukhealthcare.com.

GBUK Group Ltd.

Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0) 1757 288 587
US: 18665 289 422

www.gbukhealthcare.com

info@gbukhealthcare.com

QualRep Services BV

Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

Swiss AR Services AG

Industriestrasse 47,
6300 Zug,
Switzerland