**GBUK Healthcare**  
Part of the GBUK Group ●●●●

DASH 6™

Hub Cap

(BNS)

**INSTRUCTIONS FOR USE**

**GBUK Group Ltd.**  
Woodland House,  
Blackwood Hall Business Park,  
North Duffield, Selby,  
North Yorkshire, YO8 5DD, UK  
**EU: +44 (0) 1757 288 587**  
**US: 18665 289 422**  
[www.gbukhealthcare.com](http://www.gbukhealthcare.com)  
[info@gbukhealthcare.com](mailto:info@gbukhealthcare.com)

**EC REP** **QualRep Services BV**  
Utrechtseweg 310 B42  
6812 AR Arnhem,  
The Netherlands

**CH REP** **Swiss AR Services AG**  
Industriestrasse 47,  
6300 Zug,  
Switzerland

	<b>Page</b>
<b>EN</b> DASH 6™ Hub Cap	4
<b>FR</b> Capuchon D'embase DASH 6™	5
<b>IT</b> Cappuccio Per Mozzo DASH 6™	6
<b>ES</b> DASH 6™ Hub Cap	7
<b>DE</b> DASH 6™ Nabenerschlusskappe	8
<b>DA</b> DASH 6™ Navhætte	9
<b>NO</b> DASH 6™ Navhette	10
<b>SE</b> DASH 6™ Navlock	11
<b>NL</b> DASH 6™-Hubdop	12
<b>PL</b> Zatyczka Do Nakładki DASH 6™	13
<b>PT</b> Tampa Do Núcleo DASH 6™	14
<b>EL</b> DASH 6™ Καλυπτρα Πλημνης	15
<b>TR</b> DASH 6™ Göbek Kapağı	16
<b>CS</b> DASH 6™ – Uzávěr Objímky	17
<b>ET</b> DASH 6™ Rummukork	18
<b>HU</b> DASH 6™ Fecskendőtükrapak	19
<b>LV</b> DASH 6™ Rumbas Vāciņš	20
<b>MT</b> Tapp Tal-Buttun DASH 6™	21
<b>RO</b> DASH 6™ Capac Butuc	22
<b>SL</b> Pokrovček DASH 6™	23
<b>SK</b> Krytka Hlavice DASH 6™	24
<b>BG</b> Капачка За Главина DASH 6™	25
<b>LT</b> DASH 6™ Įvorės Gaubtelis	26

EN

## DASH 6™ HUB CAP

### Instructions For Use

Please read these instructions before use.

### INTENDED USE

The DASH 6™ Hub Cap is intended for sealing the hub of an ISO 80369-6 NRFit® compliant needle or accessory.

### INDICATION

The device is used with an ISO 80369-6 NRFit® compliant neuraxial needle or accessory and is used to seal the port when the device is not in use.

### INTENDED USER

The device is intended to be used by a trained healthcare professional.

### CONTRAINDICATIONS

The device is not intended for parenteral or enteral use.

### WARNINGS & PRECAUTIONS

- Do not use if packaging is damaged.
- Do not re-use or reprocess the device as this may lead to leakage or infection.
- The device is single use only. Reuse may result in loss of device performance and/or infection.

### DIRECTIONS

Always follow the hospital protocol when using the device.

The device is to be used in conjunction with the ISO 80369-6 NRFit® devices. Please follow the instructions for use as per the ISO 80369-6 NRFit® device provided for sealing the port.


Attach the device to an ISO 80369-6 NRFit® needle or accessory by twisting the cap onto the ISO 80369-6 NRFit® compliant device.

### STORAGE CONDITION

Keep dry and away from direct sunlight.

### DISPOSAL

Dispose of the device as per hospital protocol.

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GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

<sup>™</sup>Trademark or Registered Trademark of GBUK Group Ltd.

[www.gbukgroup.com](http://www.gbukgroup.com)

**Rx Only** Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Symbol glossary available at [www.gbukhealthcare.com](http://www.gbukhealthcare.com).

 **GBUK Group Ltd.**  
Woodland House,  
Blackwood Hall Business Park,  
North Duffield, Selby,  
North Yorkshire, YO8 5DD, UK  
**EU: +44 (0) 1757 288 587**  
**US: 18665 289 422**  
[www.gbukhealthcare.com](http://www.gbukhealthcare.com)  
[info@gbukhealthcare.com](mailto:info@gbukhealthcare.com)

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