DASH 6™

Epidural Flat

Filter (BNS)

INSTRUCTIONS FOR USE

GBUK Group Ltd.
Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0) 1757 288 587
US: 18665 289 422
www.gbukhealthcare.com
info@gbukhealthcare.com

EC REP **QualRep Services BV**
Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

CH REP **Swiss AR Services AG**
Industriestrasse 47,
6300 Zug,
Switzerland

	Page
EN DASH 6™ Epidural Flat Filter	4
FR Filtre Épidural Plat DASH 6™	5
IT Filtro Piatto Epidurale DASH 6™	6
ES Filtro Epidural Plano DASH 6™	7
DE DASH 6™ Epiduralfachfilter	8
DA DASH 6™ Epiduralt Fladt Filter	9
NO DASH 6™ Epiduralt Flatt Filter	10
SE DASH 6™ Epiduralt Platt Filter	11
NL DASH 6™ Epidural Flat Filter	12
PL Filtr Płaski Zewnątrzoponowy DASH 6™	13
PT Filtro Epidural Plano DASH 6™	14
EL Επιπεδο Επισκληριδιο Φίλτρο DASH 6™	15
TR DASH 6™ Epidural Düz Filt	16
CS DASH 6™ – Epidurální Plochý Filtr	17
ET DASH 6™ Epiduraalne Lamefilter	18
HU DASH 6™ Epidurális Lapos Szűrő	19
LV DASH 6™ Epidurālais Plakanais Filtrs	20
MT Filtru Ċatt Epidurali DASH 6™	21
RO DASH 6™ Filtru Plat Epidural	22
SL Epiduralni Ploski Filter DASH 6™	23
SK Epidurálny Plochý Filter DASH 6™	24
BG Епидурален Плосък Филтър DASH 6™	25
LT DASH 6™ Plokščiasis Epidurinis Filtras	26

EN

DASH 6™ EPIDURAL FLAT FILTER

Instructions For Use

Please read these instructions before use.

INTENDED USE

The DASH 6™ Epidural Flat Filter is intended to be used with ISO 80369-6 NRFit® compliant devices to ensure aseptic administration of neuraxial medication and anaesthetic.

INDICATION

The device is used with ISO 80369-6 NRFit® compliant neuraxial devices to maintain sterility of neuraxial medication.

INTENDED USER

The device is intended to be used by a trained healthcare professional.

CONTRAINDICATIONS

The device is not intended for parenteral or enteral use.

WARNINGS & PRECAUTIONS

Do not use if packaging is damaged.

Do not re-use or reprocess the device as this may lead to leakage.

The device is single use only. Reuse may result in loss of device performance and/or infection.

DIRECTIONS

Always follow the hospital protocol when using the device.

The device is to be used in conjunction with an ISO 80369-6 NRFit® syringe and epidural access devices/needles. Please follow the instructions for use as per the hospital procedure/instruction when using the syringe and epidural access device/needle.

The Epidural Flat Filter is connected in line between the syringe and the epidural access device/needle.

STORAGE CONDITION

Keep dry and away from direct sunlight.

DISPOSAL

Dispose of the device as per hospital protocol.

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Notes

GBUK Group award-winning specialist companies are experts in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

[™]Trademark or Registered Trademark of GBUK Group Ltd.

www.gbukgroup.com

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Symbol glossary available at www.gbukhealthcare.com.

GBUK Group Ltd.

Woodland House,
Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire, YO8 5DD, UK
EU: +44 (0) 1757 288 587
US: 18665 289 422

www.gbukhealthcare.com

info@gbukhealthcare.com

QualRep Services BV

Utrechtseweg 310 B42
6812 AR Arnhem,
The Netherlands

Swiss AR Services AG

Industriestrasse 47,
6300 Zug,
Switzerland