

SLING/HOIST COMPATIBILITY STATEMENT

Standard loop-strap sling products within the Care & Independence GLOVE SLING® and COMFORT VALUE® ranges are deemed compatible with the hoist models from the manufacturers' listed below.

In the event of any doubt or if a specific hoist does not appear in the list, please contact us to confirm compatibility *before* using the Care & Independence sling.

NOTE: Care & Independence slings are manufactured as standard with a loop strap configuration. In some instances, a clip variant may be available upon request. Please contact us for further details. In the event of having a Care & Independence clip sling variant, please ensure the sling is used with an appropriate spreader bar purposely designed for clip attachments.

MANUFACTURER	MODEL TYPES	
	Track Hoists	Mobile Hoists
Alerta		Powerlifter 150 & 175
Arjo	Maxi Sky-2, 2-plus, 1000 & 440	Maxi -500, Move, Twin, Minstrel & Tenor
Aldersly		Bariatric Mobile Hoist
BHM	Voyager, V4, Maxi	
Chiltern	Wispa 100, 200, 300, Wispa Lite	Wispa Highlift Mobile Hoist
Days Medical DMA		Mini, Maxi
Drivedevilbiss		Samssoft 150, 175. Casa Lift 180E
Ergolet	Luna	
Gulldman	Dan-Hoist/BASIC. GH2, GH3.	GL5
Handicare	Rise Atlas, basic track hoists	Carina, Eva, Vega
Handi-Move		Victor
Helping Hand	Skyframe	Swan
Homecraft		Mobile hoist
Horcher	Diana* & Lexa	<i>*Excludes GHUG Mk1</i>
Huntleigh		Porta 127, 150, 155, 175 and 225
Invacare	Robin, 800	Birdie, Alto, Compact, Reliant and Saxo
Innova	AirRise 200, 275, 350 & 500	Mobile Rise 250
Joerns/Oxford	Voyager 420, 550 & 800 Voyager 420, 550, 800, Elara & Voyager Portable	Advance, Presence, Mini, Midi & Maxi
Liftech	Gemini	
Liko	Likorail, Multirail & Ultra Twin	Golvo, Viking and Uno
Macworth	Essential. M Series	
Molift	Air, Nomad	Smart, Mover & Partner
Monarch	Track hoist	
Morecare	Swift Lifter	
Nopac	Windsor, Royal	
Prism Medical	Freeway Transactive, Prism CP200	Freeway M Series, S Series, A Series, Freeway P200
Smith & Nephew		Liftmaster
Southern Care	H130, H132, H190, HS200, EasyHoist	
Symetrik	Skyframe	
Unihoist Ltd		University Scholar, Graduate, Tutor and Master
Test Valley		Nova
Tri-Care	Freeway mobile, track and portable hoist	
Waverly Glen	Transactive, Transportable track hoist	
Wessex	Travelmaster	