

Materni

Teats

Materni's ready to use range of teats for neonatology, maternity & paediatrics. Designed for functionality.

Accessories


Teats

Materni teats grow with your baby. There are suitable sizes for every baby, for both standard and thicker milks, and across all age ranges.

Finding the perfect teat for each baby is dependent on the baby's age, latching skills and consistency of the preferred milk. Materni offers a broad range of affordable teats with different sizes and speeds, suitable for every situation.

- ✓ Safe and hygienic, individually packed and ready to use
- ✓ Zero chemicals philosophy: ETO-, Latex-, BPA-, phthalate and volatiles free.
- ✓ Bite resistant compliant (EN 14350:2004 (1&2) – (EU) n° 10-2011)
- ✓ No teat collapse possible during feeding because of reinforced design
- ✓ Ergonomic: wide teat base
- ✓ Fit for different ages: small (premature), medium (< 2.75 kg) large (> 2.75 kg)
- ✓ Wide range: 50, 80, 120 and 240 ml
- ✓ Choose a desired flow function:
 - ✓ Slow Flow (SF) – hard working – breast feed supportive – thin milk
 - ✓ Fast Flow (FF) – for thicker milk
 - ✓ Extra Fast Flow (XFF) for high viscous formula milk
- ✓ Universal fit for all feeding bottles
- ✓ Unique one-piece ring & teat design = no leak

Materni teats have unique advantages:



Small teat

Neonatology



Slow Flow

When the baby and I arrived home, the proud dad had already spent 60 euros on luxury branded teats and bottles, but none of them were satisfactory. We could not stop the dripping. Fortunately, after just a quick look on the internet, we found Materni hospital teats

“

”

Mother Karen

Medium teats

Maternity & Paediatrics



Slow Flow

Fast Flow

Extra Fast Flow

Large teats

Maternity & Paediatrics



Slow Flow

Fast Flow

Extra Fast Flow





Materni safe & easy standard product range RTU bottle and closure

Direct Codes	Product Descriptions	Screw Thread	Microbiology	Bundles/Box	Items/Box
MEPA0101	Materni teat - Large - Slow Flow	EU	CBM**	Individually packed	140
MEPA0102	Materni teat - Large - Fast Flow	EU	CBM**	Individually packed	140
MEPA0115	Materni teat - Large - Extra Fast Flow** (NEW)	EU	CBM**	Individually packed	140
MEPA0103	Materni teat - Medium - Slow Flow	EU	CBM**	Individually packed	140
MEPA0104	Materni teat - Medium - Fast Flow	EU	CBM**	Individually packed	140
MEPA0118	Materni teat - Medium - Extra Fast Flow** (NEW)	EU	CBM**	Individually packed	140
MEPA0105	Materni teat - Neonatology / Small - Slow Flow	EU	CBM**	Individually packed	140

* CBM: Clean by manufacturing, clean room ISO 6 produced and packed. Microbiologically clean.

** For thick (viscous) milk types



GBUK Group award-winning specialist companies are expert in the provision of enteral feeding, critical care and patient moving and handling products and technologies. GBUK Group companies are suppliers to every UK NHS Trust and to a rapidly expanding export customer base.

www.gbukgroup.com



www.gbukenteral.com



GBUK Enteral Ltd.
Woodland House
Blackwood Hall Business Park
North Duffield, Selby
North Yorkshire YO8 5DD

T +44 (0) 1757 282 945
E info@gbukenteral.com



Nursing innovations
Clean by manufacturing

www.maternimedical.eu