

## Patient Handling Devices

### GBUK Banana

Woodland House, Blackwood Hall Business Park  
North Duffield, Selby, North Yorkshire YO8 5DD

T +44 (0) 1757 288 587  
W [www.gbukbanana.com](http://www.gbukbanana.com)  
E [info@gbukbanana.com](mailto:info@gbukbanana.com)

 @GBUK\_Group  
 @GBUKBanana  
 GBUK Group  
 GBUK Group



## GBUK Banana Clip Sling Compatibility Statement

### Applicable product codes:

B-SSPUC-XS  
B-SSPUC-S  
B-SSPUC-M  
B-SSPUC-L  
B-SSPUC-XL

### Quality

GBUK Banana Clip Slings have been designed and tested in compliance with ISO 10535:2006 “Hoists for the Transfer of Disabled Persons- Requirements and Test Methods”.

The GBUK Banana Clip Slings are CE marked and are a Class I Medical Device.

FDA Device Listings Number D422414.

GS1 compliant artwork on packaging.

### Compatibility

Hoists and slings are widely interchangeable if both are risk assessed prior to use and confirmed to be compatible.

GBUK Banana have carried out compatibility tests on the following limited brands of hoists currently available to us at this time and found that these are compatible with the GBUK Banana Single Patient Use Clip Slings.

Oxford Joerns Presence  
Oxford Joerns Advance  
Oxford Adjustable Power Cradle  
Arjo Maxi Move

If the hoist you use is not listed above, we would be pleased to complete a compatibility assessment and statement for the brand you are using. Please contact GBUK Banana via the email address below and we will be happy to arrange this.

GBUK Banana always recommend you conduct a risk assessment prior to use, and LOLER regulations require you to carry out a visual inspection of the sling prior to each use.

### Contact Information

GBUK Banana  
Woodland House,  
Blackwood Hall Business Park,  
North Duffield,  
Selby,  
North Yorkshire  
YO8 5DD  
Telephone: +44 (0) 1757 288 587  
Email: [info@gbukbanana.com](mailto:info@gbukbanana.com)  
Website: [www.gbukbanana.com](http://www.gbukbanana.com)

