

Grip-Lok[®] 2100NGH Neonatal Nasal Gastric Tube Securement

- ✓ Hypoallergenic, breathable and latex-free skin contact adhesive
- ✓ Hydrocolloid adhesive is water-soluble pectin based, designed to be very gentle on fragile skin
- ✓ Soft and flexible fabric design improves patient comfort
- ✓ Easy to apply, open, inspect and adjust
- ✓ Provides superior securement for both vertical and horizontal lifting accidental line pulls
- ✓ Can easily be placed and adjusted with gloved hands
- ✓ Sterile packed

Grip-Lok[®] Neonatal Nasal Gastric offers firm securement of feeding tubes including pure silicone materials and will fit a wide range of tubes.



Note: The Grip-Lok 2100NGH may be cut in half prior to removal of liners to reduce the size



Reference	NPC	Description	Tubing Range	Device Length	Units/ Box
2100NGH	ELW931	Neonatal Hydrocolloid Nasal Gastric Tube Securement	1.0mm - 2.0mm	4.2cm	50 / 1

