


# Grip-Lok<sup>®</sup> 2100ANG Nasal Gastric Securement

- ✓ Hypoallergenic, breathable and latex-free skin contact adhesive
- ✓ Provides superior securement for both vertical and horizontal lifting accidental line pulls
- ✓ Split 3-section lid for enhanced tube securement
- ✓ Side notch to assist with correct alignment against nostril
- ✓ Soft and flexible fabric design
- ✓ Easy to apply, open, inspect and adjust
- ✓ Can easily be placed and adjusted with gloved hands
- ✓ Sterile packed

**Grip-Lok<sup>®</sup> Nasal Gastric offers firm securement of feeding tubes and will fit a wide range of tube sizes**



Reference	NPC	Description	Tubing Range	Device Length	Units/ Box
2100ANG	ELW705	Small Tube - Nasal Gastric tube securement	3.0mm - 5.0mm	3.8cm	100 / 1