

Suction Catheter Size Guide


Updated hub colour change in line with the ISO 8836:2019 Suction Catheters for use in the Respiratory tract standard

Always check product label if in doubt of size


| Size (Fr) | 5 | 6 | 7 | 8 | 10 | 12 | 14 | 16 |
|-----------|------|-------------|-------|------------|-------|-------|-------|--------|
| Colour | Grey | Light Green | Ivory | Light Blue | Black | White | Green | Orange |




Please be aware 7Fr has changed from purple to ivory in January 2022




| | | |
|------------|--------------------------------|---------------|
| 5Fr | 100 pcs / box Product Code: | NHSSC: |
| 35cm | TT01-05-035 | FSQ587 |
| 50cm | TT01-05-050 | FSQ993 |



| | | |
|------------|--------------------------------|---------------|
| 6Fr | 100 pcs / box Product Code: | NHSSC: |
| 35cm | TT01-06-035 | FSQ588 |
| 50cm | TT01-06-050 | FSQ568 |




| | | |
|------------|--------------------------------|---------------|
| 7Fr | 100 pcs / box Product Code: | NHSSC: |
| 35cm | TT01-07-035 | FSQ589 |
| 50cm | TT01-07-050 | FSQ570 |




| | | |
|------------|--------------------------------|---------------|
| 8Fr | 100 pcs / box Product Code: | NHSSC: |
| 35cm | TT01-08-035 | FSQ591 |
| 50cm | TT01-08-050 | FSQ572 |




| | | |
|-------------|--------------------------------|---------------|
| 10Fr | 100 pcs / box Product Code: | NHSSC: |
| 35cm | TT01-10-035 | FSQ593 |
| 50cm | TT01-10-050 | FSQ588 |
| 60cm | TT01-10-060 | FSQ580 |



| | | |
|-------------|--------------------------------|---------------|
| 12Fr | 100 pcs / box Product Code: | NHSSC: |
| 50cm | TT01-12-050 | FSQ570 |
| 60cm | TT01-12-060 | FSQ575 |



| | | |
|-------------|--------------------------------|---------------|
| 14Fr | 100 pcs / box Product Code: | NHSSC: |
| 60cm | TT01-14-060 | FSQ578 |



| | | |
|-------------|--------------------------------|----------------|
| 16Fr | 100 pcs / box Product Code: | NHSSC: |
| 60cm | TT01-16-060 | FSQ1787 |