

Neuraxial Devices

NRFit™ Syringes, Needles & Accessories



Improving patient safety in neuraxial connections.

ISO 80369-6 NRFit™ neuraxial connectors

ISO 80369-6:2016

is the international standard that specifies requirements for small-bore connectors intended to be used for connections in neuraxial applications.

Neuraxial Applications

involve the use of medical devices intended to administer medications to neuraxial sites, wound infiltration, anaesthesia delivery, and other regional anaesthesia procedures or to monitor or remove cerebro-spinal fluid for therapeutic or diagnostic purposes.

NRFit™

is the dedicated connector for neuraxial and regional block devices and is not compatible with luer connectors.

Dash 6™

is the exclusive brand name for the NRFit™ connectors from GBUK Healthcare.



Lock & Slip Syringes:

Product Code	Description	Size	Quantity
15011010	Lock Syringe	1ml	100 pcs/box
15011030	Lock Syringe	3ml	100 pcs/box
15011050	Lock Syringe	5ml	100 pcs/box
15011100	Lock Syringe	10ml	100 pcs/box
15011200	Lock Syringe	20ml	50 pcs/box
15011600	Lock Syringe	60ml	25 pcs/box



Product Code	Description	Size	Quantity
15021010	Slip Syringe	1ml	100 pcs/box
15021030	Slip Syringe	3ml	100 pcs/box
15021050	Slip Syringe	5ml	100 pcs/box
15021100	Slip Syringe	10ml	100 pcs/box
15021200	Slip Syringe	20ml	50 pcs/box



Accessories:

Product Code	Description	Quantity
15241001	3-Way Stopcock	50 pcs/box
15111001	Needle Hub Cap	100 pcs/box
15041001	Syringe Cap	100 pcs/box
15171001	Syringe to Syringe Adaptor	100 pcs/box
15231022	0.22µ Disc Filter	100 pcs/box
15311022	0.22µ Epidural Flat Filter	100 pcs/box
15281100	Drawing Up Filter Straw 5µ 100mm	100 pcs/box
15491050	Blunt Drawing Up Needle without Filter 18G x 50mm	100 pcs/box
15291050	Blunt Drawing Up Needle with 5µ Filter 18G x 50mm	100 pcs/box
15461010	Extension Set with Clamp 10cm	25 pcs/box
15461020	Extension Set with Clamp 20cm	25 pcs/box
15461125	Extension Set with Clamp 125cm	25 pcs/box
15351225	Butterfly Needle 25G	100 pcs/box
15351227	Butterfly Needle 27G	100 pcs/box
15361221	Butterfly Needle 21G with one-way valve	100 pcs/box

