



INSTRUCTION FOR USE:

GBUK Safety Syringe with Fixed Needle

A.




B. Hold here while removing cap



C. After use, push cover over needle
until a click is heard



D. Hold barrel, twist safety cover to lock,
until a click is heard



GBUK Group Ltd.

Woodland House, Blackwood Hall Business Park
North Duffield, Selby, North Yorkshire, YO8 5DD, United Kingdom

Email: info@gbukhealthcare.com
Main: +44 (0)1757 288 587

www.gbukgroup.com/safety-syringe-with-fixed-needle/



GBUKHealthcare
Part of the GBUK Group