



Comfort
VALUE™

For stores and catalogues

COMFORT VALUE™ RECLINE HIGHBACK (CVR-HB)

PART OF THE IN-CHAIR SLING FAMILY

A long-seat sling with fixed pommel and integrated head support making the COMFORT VALUE™ Recline Highback a great choice as a stocked product.

FEATURES & BENEFITS

- Designed for in-situ use
- Full body support
- Fixed pommel and padded roll adding comfort and removing stress from between legs and knees
- Outside leg pickup
- Structured sling shape via darting but without seams for user comfort
- Convenient pockets for tucking in pommel tapes and leg loop straps when not in use
- Colour-coded binding for quick size recognition
- Suitable for general transfer as well as hoisting user in supine position
- This Highback variant has an integrated head support that can be folded down and used as a neck support.
- Suitable for 95% of users - including amputees



CVR-HB



CVR-HB

SLING SHAPE	In-Chair
PRIMARY USE	
TRANSFER	
SIZE	
SWL (KG)	
ATTACHMENT TYPE	3-Loop strap
CLIP OPTION	✗
MODESTY COVER	Pommel
PICK UP	Outside leg
POSITIONING AID	✗
HEAD SUPPORT	✓
APERTURE	✗
AMPUTEE	✓
MATERIALS	3D-Spacer
PADDING	✗
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✗
WARRANTY	3 year
COMPATIBILITY	Majority of hoists
OTHER INFO	Pockets for straps when not in use Non-assessed product for stores & catalogues



Call 033 012 80808 to discuss your requirements with a Business Development Manager