



Comfort
VALUE™

For stores and catalogues

COMFORT VALUE™ RECLINE (CVR)

PART OF THE IN-CHAIR SLING FAMILY

A long-seat sling with fixed pommel that makes the COMFORT VALUE™ Recline slings a great stocked product.

FEATURES & BENEFITS

- Designed for in-situ use
- Full body support
- Fixed pommel and padded roll adding comfort and removing stress from between legs and knees
- Outside leg pickup
- Structured sling shape via darting but without seams for user comfort
- Convenient pockets for tucking in pommel tapes and leg loop straps when not in use
- Colour-coded binding for quick size recognition
- Suitable for general transfer as well as hoisting user in supine position
- Suitable for 95% of users - including amputees

Also available as Highback variant



CVR



CVR

SLING SHAPE	In-Chair
PRIMARY USE	
TRANSFER	
SIZE	XXS-XL
SWL (KG)	230
ATTACHMENT TYPE	3-Loop strap
CLIP OPTION	✗
MODESTY COVER	Pommel
PICK UP	Outside leg
POSITIONING AID	✗
HEAD SUPPORT	✗
APERTURE	✗
AMPUTEE	✓
MATERIALS	3D-Spacer
PADDING	✗
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✗
WARRANTY	3 year
COMPATIBILITY	Majority of hoists
OTHER INFO	Pockets for straps when not in use Non-assessed product for stores & catalogues



Call 033 012 80808 to discuss your requirements with a Business Development Manager