



Comfort
VALUE™

For stores and catalogues

COMFORT VALUE™ STAND ASSIST (CVSA)
THE SPECIALIST SLING FAMILY

The COMFORT VALUE™ Stand Assist is a useful aid to help move from a sit-to-stand position and is great for single-handed care.

FEATURES & BENEFITS

- Easy to fit
- Underarm padding for user comfort
- Non-slip inner back panel
- Hook and Loop fastening chest band with integral clunk-click buckle for user safety and security



CVSA

SIT-TO-STAND

SLING SHAPE	Specialist
PRIMARY USE	
TRANSFER	
SIZE	
SWL (KG)	
ATTACHMENT TYPE	3-Loop strap
CLIP OPTION	✗
MODESTY COVER	N/A
PICK UP	Underarm
POSITIONING AID	✗
HEAD SUPPORT	N/A
APERTURE	✓
AMPUTEE	N/A
MATERIALS	Soft-knit Polyester + non-slip back panel
PADDING	Underarm
STANDARD COLOUR	
ACCESSORIES	✗
CUSTOMISABLE	✗
WARRANTY	3 years
COMPATIBILITY	Majority of hoists
OTHER INFO	Hook & Loop chest band + clunk-click buckle Non-assessed product for stores & catalogues



Call 033 012 80808 to discuss your requirements with a Business Development Manager