



Comfort
VALUE™

For stores and catalogues

COMFORT VALUE™ TRANSPORT (CVTR)

THE SPECIALIST SLING FAMILY

NEW

The COMFORT VALUE™ Transport aids in the transfer of users from a sit-to-stand position whilst remaining in a semi-seat position. Ideal for times when a full sling transfer is not required.

FEATURES & BENEFITS

- Easy to fit
- Extra padding to legs and underarms for user comfort
- Non-slip inner back panel
- Clunk-click buckle for user safety and security



CVTR



CVTR

SIT-TO-STAND

SLING SHAPE	Specialist
PRIMARY USE	
TRANSFER	
SIZE	XS-XL
SWL (KG)	230
ATTACHMENT TYPE	4-Loop strap
CLIP OPTION	✗
MODESTY COVER	N/A
PICK UP	Inside leg
POSITIONING AID	✗
HEAD SUPPORT	N/A
APERTURE	✓
AMPUTEE	✗
MATERIALS	Soft-knit Polyester + non-slip back panel
PADDING	Underarm and leg sections
STANDARD COLOUR	
ACCESSORIES	✗
CUSTOMISABLE	✗
WARRANTY	3 years
COMPATIBILITY	Majority of standing aids
OTHER INFO	Non-slip back panel + clunk-click buckle Non-assessed product for stores & catalogues



Call 033 012 80808 to discuss your requirements with a Business Development Manager