



Comfort
VALUE™

For stores and catalogues

COMFORT VALUE™ UNIVERSAL POLY SIX (CVUPS)
PART OF THE UNIVERSAL SLING FAMILY

The ‘Poly Six’ variants in the COMFORT VALUE™ Universal sling family have two additional hip straps to help prevent the user abducting during hoisting.

FEATURES & BENEFITS

- Full body support with tapered leg sections
- Features double-sided Parasilk™ slide sheet material on leg pieces for ease of fitting
- Binding-free outside leg sections to reduce pressure behind the knee and reduce shear and drag during fitting
- Soft and comfortable being fashioned from a one-cut soft-knit Polyester fabric
- Secure integrated 3-loop straps at shoulder and leg
- Additional 2-loop hip strap reduces the risk of abduction of the hips when hoisting
- Double-stitched colour-contrast seams at vulnerable areas for added durability
- Easy to apply with straight leg and inside pick-up
- Aperture for user convenience when toileting or for access
- Can be hammocked for amputee use without need for specialist sling



CVUPS



CVUPS

GENERAL TRANSFER

SLING SHAPE	Universal
PRIMARY USE	GENERAL TRANSFER
TRANSFER	
SIZE	S-XL
SWL (KG)	230
ATTACHMENT TYPE	3-Loop strap + hip strap
CLIP OPTION	✗
MODESTY COVER	Loop
PICK UP	Inside leg
POSITIONING AID	✗
HEAD SUPPORT	✗
APERTURE	✓
AMPUTEE	✓
MATERIALS	Soft-knit Polyester + doubled-sided Parasilk™ (legs)
PADDING	Legs
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✗
WARRANTY	3 year
COMPATIBILITY	Majority of hoists
OTHER INFO	Non-assessed product for stores & catalogues



Call 033 012 80808 to discuss your requirements with a Business Development Manager