



Comfort
VALUE™

For stores and catalogues

COMFORT VALUE™ UNIVERSAL POLY (CVUP)
PART OF THE UNIVERSAL SLING FAMILY

The Polyester variations of the COMFORT VALUE™ Universal sling style offer extra user comfort with binding-free leg sections that aid ease of fitting too.

FEATURES & BENEFITS

- Full body support with tapered leg sections
- Parasilk™ slide sheet material to leg section for ease of fitting
- Soft and comfortable being fashioned from a one-cut soft-knit Polyester fabric
- Secure integrated strap fixings
- Double-stitched colour-contrast seams at vulnerable areas for added durability
- Easy to apply with straight leg and inside pick-up
- Aperture for user convenience when toileting or for access
- Can be hammocked



CVUP



CVUP

GENERAL TRANSFER

SLING SHAPE	Universal
PRIMARY USE	GENERAL TRANSFER
TRANSFER	
SIZE	S-XL
SWL (KG)	230
ATTACHMENT TYPE	3-Loop strap
CLIP OPTION	✗
MODESTY COVER	Loop
PICK UP	Inside leg
POSITIONING AID	✗
HEAD SUPPORT	✗
APERTURE	✓
AMPUTEE	✓
MATERIALS	Soft-knit Polyester + double-sided Parasilk™ (legs)
PADDING	Legs
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✗
WARRANTY	3 year
COMPATIBILITY	Majority of hoists
OTHER INFO	Non-assessed product for stores & catalogues



Call 033 012 80808 to discuss your requirements with a Business Development Manager