



GLOVE™ CONVENIENCE II (GC2)
PART OF THE TOILETING SLING FAMILY

Highly supportive, super secure and with a unique ‘anti-jack-knife’ design feature, our GLOVE™ Convenience II is a fantastic toileting and dressing solution - even for users previously considered unsuitable for toileting/access slings.

FEATURES & BENEFITS

- Anti-jack-knife design prevents twisting
- Adjustable dual buckle chest band
- Soft-knit Polyester with ‘Comfort’ lining
- Large aperture for ease of access
- Deep lumbar support
- Deep underarm padding for increased comfort
- Double-sided high-grade slidesheet fabric on leg and vertebrae sections
 - Minimises shear and drag
 - Reduces tissue viability and pressure injury risks



GC2



SLING SHAPE	Access
PRIMARY USE	
TRANSFER	
SIZE	3-9
SWL (KG)	3-9: 200
ATTACHMENT TYPE	3-Loop strap
CLIP OPTION	✓
MODESTY COVER	✗
PICK UP	Inside leg
POSITIONING AID	Positioning Handle (1)
HEAD SUPPORT	✗
APERTURE	✓
AMPUTEE	✗
MATERIALS	Soft-knit Polyester, slidesheet (legs) + Terryloth ‘Comfort’ lining (back panel)
PADDING	Legs, underarms and back
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	Majority of hoists



TOILETING / ACCESS

Call 033 012 80808
to book a free assessment