



GLOVE™ CONVENIENCE MAX (GCM)
PART OF THE TOILETING SLING FAMILY

The GLOVE™ Convenience Max is an elegant advanced toileting sling that encompasses a self-adjusting corset-style body with patented binding-free padded leg section offering users maximum comfort and maximum security during use. Ideal for users with minimal upper body strength!

FEATURES & BENEFITS

- Self-adjusting corset-style body design that hugs as it lifts for maximum support and comfort
- Adjustable shoulder straps
- Adjustable triple clunk-click buckle chest band
- Large aperture for maximum access and hygiene
- Padded lumbar support with luxurious one-way velvet velour lining
- Deep underarm padding for increased comfort
- Padded with 3D-Spacer fabric for comfort, durability and quicker drying
- High back support enables lifting from bed and floor as well as seated position
- Patented binding-free double-sided slidesheet padded leg sections aids fitting and reduces skin sheer and tear
- Positioning Pockets for ease of fitting, greater user dignity during hoisting and better safe patient handling practice



GCM



GCM

SLING SHAPE	Access
PRIMARY USE	
TRANSFER	
SIZE	4-9
SWL (KG)	4-9: 230
ATTACHMENT TYPE	Adjustable shoulder strap and 3-Loop leg strap
CLIP OPTION	✗
MODESTY COVER	✗
PICK UP	Inside leg
POSITIONING AID	Positioning Pocket (2)
HEAD SUPPORT	✗
APERTURE	✓
AMPUTE	✗
MATERIALS	Soft-knit Polyester, slidesheet (legs) One-way Velvet Velour
PADDING	3D-Spacer to lumbar, underarm and legs
STANDARD COLOUR	
ACCESSORIES	✗
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	Majority of hoists
OTHER INFO	Triple clunk-click buckle chestband



TOILETING / ACCESS

Call 033 012 80808
to book a free assessment