



GLOVE™ CONVENIENCE PLUS (GCP) PART OF THE TOILETING SLING FAMILY

The GLOVE™ Convenience Plus sling provides a toileting solution with great access for the carer and security for the user.

FEATURES & BENEFITS

- Anti-jack-knife design prevents twisting
- Wrap-over chest band with integrated adjustable buckle
- Soft-knit Polyester with 'Comfort' lining
- Large aperture for ease of access
- Double-sided high-grade slidesheet fabric on leg and vertebrae sections
 - Minimises shear and drag
 - Reduces tissue viability and pressure injury risks
 - Soft to the touch, extra slippery and quiet application



GCP



GCP

SLING SHAPE	Access
PRIMARY USE	
TRANSFER	
SIZE	3-9
SWL (KG)	3-9: 200
ATTACHMENT TYPE	Adjustable Extensor strap
CLIP OPTION	✓
MODESTY COVER	✗
PICK UP	Outside leg
POSITIONING AID	Positioning Handle (3)
HEAD SUPPORT	✓
APERTURE	✓
AMPUTEE	✗
MATERIALS	Soft-knit Polyester + slidesheet (legs) + Terrycloth 'Comfort' lining to body
PADDING	Lumber
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	Majority of hoists



**Call 033 012 80808
to book a free assessment**