



GLOVE™ DELUXE BATHING (GXB)
PART OF THE DELUXE SLING FAMILY

The GLOVE™ Deluxe Bathing is comfortable, supportive and quick draining making it a great choice for those with bathing support needs.

FEATURES & BENEFITS

- Soft-knit Polyester for improved tissue viability
- 3D-Polyester in leg sections
 - faster drainage
 - added user comfort
- Modesty loop helps prevent user slipping in wet environment



GXB



GXB

SLING SHAPE	Deluxe
PRIMARY USE	 
TRANSFER	   
SIZE	 2-12
SWL (KG)	2-9:  230 10-12:  350
ATTACHMENT TYPE	3-loop strap
CLIP OPTION	✓
MODESTY COVER	Loop
PICK UP	Outside leg
POSITIONING AID	Positioning Handle (2)
HEAD SUPPORT	✗
APERTURE	✓
AMPUTEE	✓
MATERIALS	Soft-knit Polyester
PADDING	Legs
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	Majority of hoists



BATHING

Call 033 012 80808
to book a free assessment