



GLOVE™ DELUXE TUX PLUS (GXTP)

PART OF THE DELUXE SLING FAMILY

A superb all-rounder and one of our most popular general transfer slings.

Made from a soft-knit Polyester with additional padded support, the GLOVE™ Deluxe Tux Plus provides superb comfort for users.

FEATURES & BENEFITS

- Easy to fit and remove
- Double-sided high-grade slidesheet fabric on leg and vertebrae sections
 - Minimises shear and drag
 - Reduces tissue viability and pressure injury risks
 - Soft to the touch, extra slippery and quiet application
- Great for heavier individuals
- Padded leg and vertebrae support
- Extensor strap fitted as standard
 - Aids user stability during hoist through unique "suspension" facility that counterpoises extensor spasm
 - Simplifies attachment to the spreader bar
- Integral head support
- Slightly longer leg section to support larger bodies
- Suitable for almost all conditions and disabilities



GXTP



GXTP

GENERAL TRANSFER

SLING SHAPE	Deluxe
PRIMARY USE	GENERAL TRANSFER
TRANSFER	
SIZE	
SWL (KG)	2-9: 10-12:
ATTACHMENT TYPE	Extensor strap
CLIP OPTION	✓
MODESTY COVER	Loop
PICK UP	Outside leg
POSITIONING AID	Positioning Handle (3)
HEAD SUPPORT	✓
APERTURE	✓
AMPUTEE	✓
MATERIALS	Soft-knit Polyester + slidesheet (legs)
PADDING	Legs, vertebrae and head support
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	Majority of hoists



Call 033 012 80808
to book a free assessment