



GLOVE™ DIAMOND DELUXE FOXY (GDDF)

PART OF THE DIAMOND DELUXE SLING FAMILY

Introducing the GLOVE™ Diamond Deluxe Foxy (GDDF).

With wider leg and hip sections, his sling is ideally suited for pear-shaped users with broader hips. It also benefits from all the design and performance elements associated with the Diamond Family.

FEATURES & BENEFITS

- **Patented** diamond-shaped aperture reduces pressure on the sacrum and coccyx when hoisted and provides access for toileting
- **Patented** binding-free leg sections for additional user comfort
- Double-sided high-grade slidesheet fabric on leg and vertebrae sections
 - Minimises shear and drag
 - Reduces tissue viability and pressure injury risks
 - Soft to the touch, extra slippery and quiet application
- 3-Loop straps for adjustable angle of recline
- Modesty loop for increased user security and discretion
- Padded with 3D-Spacer fabric for comfort, durability and quicker drying
- Positioning Pockets for ease of fitting, greater user dignity during hoisting and better safe patient handling practice
- Outside leg pick-up for greater control and thigh support
- Suitable for almost all conditions and disabilities
- Suitable for Amputees (above-knee: single or bilateral)



GDDF

SLING SHAPE	Deluxe
PRIMARY USE	GENERAL TRANSFER
TRANSFER	
SIZE	4-9
SWL (KG)	4-9: 230
ATTACHMENT TYPE	3-Loop strap
CLIP OPTION	✓
MODESTY COVER	Loop
PICK UP	Outside leg
POSITIONING AID	Positioning Pockets (2)
HEAD SUPPORT	✗
APERTURE	✓
AMPUTEE	✓
MATERIALS	Soft-knit Polyester + slidesheet
PADDING	Legs & vertebrae
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	Majority of hoists



GENERAL TRANSFER

Call 033 012 80808
to book a free assessment