



GLOVE™ DIAMOND DELUXE PLUS (GDDP)

PART OF THE DIAMOND DELUXE SLING FAMILY

The new GLOVE™ Diamond sling family is unrivalled for comfort and performance and when you need a little extra leg length and an integral head support, Plus variants are the ideal solution!

FEATURES & BENEFITS

- **Patented** diamond-shaped aperture reduces pressure on the sacrum and coccyx when hoisted and provides access for toileting
- **Patented** binding-free leg sections for additional user comfort
- Double-sided high-grade slidesheet fabric on leg and vertebrae sections
 - Minimises shear and drag
 - Reduces tissue viability and pressure injury risks
 - Soft to the touch, extra slippery and quiet application
- Integral head support
- Slightly longer leg section to support larger bodies
- Modesty loop for increased user security and discretion
- Padded with 3D-Spacer fabric for comfort, durability and quicker drying
- Extensor strap fitted as standard i) Aids user stability during hoist through unique “suspension” facility that counterpoises extensor spasm ii) Simplifies attachment to the spreader bar
- Positioning Pockets for ease of fitting, greater user dignity during hoisting and better safe patient handling practice
- Outside leg pick-up for greater control and thigh support
- Suitable for almost all conditions and disabilities



GDDP



GDDP

SLING SHAPE	Deluxe
PRIMARY USE	GENERAL TRANSFER
TRANSFER	
SIZE	
SWL (KG)	1-9: 10-12:
ATTACHMENT TYPE	Extensor strap
CLIP OPTION	✓
MODESTY COVER	Loop
PICK UP	Outside leg
POSITIONING AID	Positioning Pockets (2)
HEAD SUPPORT	✓
APERTURE	✓
AMPUTEE	✓
MATERIALS	Soft-knit Polyester + slidesheet
PADDING	Legs, vertebrae and head support
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	Majority of hoists
OTHER	A 3-Loop shoulder strap version is also available (GDDP/B)



GENERAL TRANSFER

Call 033 012 80808
to book a free assessment