



GLOVE™ EASY TURN (GTE)
PART OF THE SPECIALIST SLING FAMILY

The GLOVE™ Easy Turn as the name suggests, is the perfect sling when individuals require turning in bed.

Simple yet highly effective.

FEATURES & BENEFITS

- Double-sided high-grade slidesheet fabric on leg and vertebrae sections
 - Minimises shear and drag
 - Reduces tissue viability and pressure injury risks
- - Soft to the touch, extra slippy and quiet application
- Durable
- High performance
- Dual Hook & Loop straps for quick and easy attachment to bedframe



Easy, stress-free turning solution



GTE



SLING SHAPE	Specialist
PRIMARY USE	 For turning and personal care
TRANSFER	N/A
SIZE	 One size
SWL (KG)	 200
ATTACHMENT TYPE	Loop strap
CLIP OPTION	✓
MODESTY COVER	N/A
PICK UP	N/A
POSITIONING AID	N/A
HEAD SUPPORT	N/A
APERTURE	✓
AMPUTEE	✓
MATERIALS	Soft-knit Polyester + slidesheet
PADDING	✗
STANDARD COLOUR	
ACCESSORIES	N/A
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	Majority of hoists



MANUAL HANDLING

Call 033 012 80808
to book a free assessment