



GLOVE™ EURO MANUAL (GEURM) PART OF THE SPECIALIST SLING FAMILY

When fast fitting and ease of use is essential, the GLOVE™ Euro Manual is the ideal sling for emergency transfers without a hoist.

Can also be used in non-emergency situations where a hoist is not available.

FEATURES & BENEFITS

- Fast easy fit and removal
- No hoist required
- Double-sided high-grade slidesheet fabric on leg and vertebrae sections
 - Minimises shear and drag
 - Reduces tissue viability and pressure injury risks
- - Soft to the touch, extra slippy and quiet application
- Robust handles at leg and shoulder



No hoist required!



Ideal for emergency evacuations



GEURM



GEURM

SLING SHAPE	Universal
PRIMARY USE	GENERAL TRANSFER WITHOUT HOIST
TRANSFER	
SIZE	2-6
SWL (KG)	2-6: 200
ATTACHMENT TYPE	Handles for manual lift
CLIP OPTION	✗
MODESTY COVER	✗
PICK UP	Inside leg
POSITIONING AID	✗
HEAD SUPPORT	✓
APERTURE	✓
AMPUTEE	✓
MATERIALS	Soft-knit Polyester + slidesheet (legs)
PADDING	Legs
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	N/A



MANUAL HANDLING

**Call 033 012 80808
to book a free assessment**