



GLOVE™ EXTENSOR (GE)
PART OF THE DELUXE SLING FAMILY

The GLOVE™ Extensor incorporates a 'suspension' feature to counterpoise user extensor spasm.

This unique design innovation enables the sling to move with the spasm and return to its original position afterwards.

FEATURES & BENEFITS

- Easy to fit and remove
- Double-sided high-grade slidesheet fabric on leg and vertebrae sections
 - Minimises shear and drag
 - Reduces tissue viability and pressure injury risks
 - Soft to the touch, extra slippery and quiet application
- Great for individuals with extensor spasm and excessive movements



GE



GE

GENERAL TRANSFER

SLING SHAPE	Deluxe
PRIMARY USE	GENERAL TRANSFER
TRANSFER	
SIZE	
SWL (KG)	4-9:
ATTACHMENT TYPE	Extensor strap
CLIP OPTION	✓
MODESTY COVER	Loop
PICK UP	Outside leg
POSITIONING AID	Positioning Handle (3)
HEAD SUPPORT	✓
APERTURE	✓
AMPUTEE	✓
MATERIALS	Soft-knit Polyester + slidesheet
PADDING	Legs, vertebrae, head support
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	Majority of hoists



Call 033 012 80808
to book a free assessment