



GLOVE™ FOOT (GFOOT)
PART OF THE SPECIALIST SLING FAMILY

The GLOVE™ Foot is a great 'add-on' or stand alone additional leg support during fitting of sling or when applying dressings.

Easily slipped over heel with simple strap attachment direct to hoist spreader bar, the GLOVE™ Foot provides improved comfort to the thigh area and reduces excessive load for users in a leg cast. It can also help reduce back injury risk to carers during routine care giving or whilst aiding into/out of bed as well as being useful during therapy treatments.

FEATURES & BENEFITS

- Reduce leg load for improved comfort
- Helps reduce carer back injury risk
- Highly supportive
- Easy to attach and remove
- Adjustable strap for multi-positioning of leg
- Offered in three sizes
- Compatible with all sling types



GFOOT

LEG SUPPORT

SLING STYLE	Specialist
PRIMARY USE	Leg support
TRANSFER	N/A
SIZE	 Paediatric / Adult / Bariatric
SWL (KG)	 75
ATTACHMENT TYPE	Adjustable strap
CLIP OPTION	✗
MODESTY COVER	N/A
PICK UP	Underfoot
POSITIONING AID	N/A
HEAD SUPPORT	N/A
APERTURE	N/A
AMPUTEE	✗
MATERIALS	Soft-knit Polyester
PADDING	✓
STANDARD COLOUR	
ACCESSORIES	✗
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	All GLOVE SLING™ variants
OTHER INFO	Supplied as single unit

Call 033 012 80808
to book a free assessment

PRODSPEC/GFOOT-EN01