




GLOVE™ IN-CHAIR LITE (GIL)
PART OF THE IN-CHAIR HAMMOCK SLING FAMILY

GLOVE™ In-Chair sling styles are great for users who need to remain seated for longer periods than normal.

The GLOVE™ In-Chair LITE provides users with a soft and temperature comfortable experience.

FEATURES & BENEFITS

- High air-permeability helps reduce sweating and associated risks
- Unique shape provides full body support
- Added waist straps to help prevent abduction of hips during hoisting
- Strap pockets to eliminate snagging risk



Pressure mapping evidence available



GIL



GIL

IN-CHAIR

SLING SHAPE	In-Chair
PRIMARY USE	
TRANSFER	
SIZE	
SWL (KG)	2-9: 10-12:
ATTACHMENT TYPE	2-Loop shoulder & waist strap, 3-Loop leg strap
CLIP OPTION	✓
MODESTY COVER	Fitted pommel
PICK UP	Outside leg
POSITIONING AID	✗
HEAD SUPPORT	✗
APERTURE	✓
AMPUTEE	✓
MATERIALS	Soft-knit Polyester
PADDING	✗
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	Majority of hoists
OTHER INFO	Strap pockets



Call 033 012 80808
to book a free assessment