




GLOVE™ IN-CHAIR (GI)
PART OF THE IN-CHAIR HAMMOCK SLING FAMILY

GLOVE™ In-Chair sling styles are great for users who need to remain seated for longer periods than normal.

The GLOVE™ In-Chair offers users practical comfort.

FEATURES & BENEFITS

- High air-permeability helps reduce sweating and associated risks
- Unique shape provides full body support
- Added waist straps to help prevent abduction of hips during hoisting
- Strap pockets to eliminate snagging risk
- 'Comfort' terrycloth lining for extra softness



Pressure mapping
evidence available



GI



GI

IN-CHAIR

SLING SHAPE	In-Chair
PRIMARY USE	
TRANSFER	
SIZE	2-12
SWL (KG)	2-9: 200 10-12: 350
ATTACHMENT TYPE	2-Loop shoulder & waist strap, 3-Loop leg strap
CLIP OPTION	✓
MODESTY COVER	Fittedommel
PICK UP	Outside leg
POSITIONING AID	✗
HEAD SUPPORT	✗
APERTURE	✓
AMPUTEE	✓
MATERIALS	Soft-knit Polyester + Terrycloth 'Comfort' lining
PADDING	Legs
STANDARD COLOUR	
ACCESSORIES	✓
CUSTOMISABLE	✓
WARRANTY	Lifetime
COMPATIBILITY	Majority of hoists
OTHER INFO	Strap pockets



Call 033 012 80808
to book a free assessment