



Enhancing patient care

Endourology Brochure

Endourology portfolio

History of GBUK	4
------------------------	----------

Access products

Guidewire	7
Ureteral Catheter	12
Dual Lumen Catheter	16
Initial Puncture Needle	17
Chiba Needle	18

Dilation products

Urethral Dilator Set	19
S Shape Dilator	20
Ureteral Balloon Dilator	21
Ureteral Dilator	22
Ureteral Access Sheath	23
Fascial Dilator	25
Amplatz Dilator	26
Nephrostomy Dilator	27

Stone managment

G-Paw	28
Tipless Stone Basket	29
Tipped Stone Basket	34
Perk Basket	37
Stone Grasper	38
Stone Stop	39
URS Forceps	40

Endourology portfolio

Drainage products

Ureteral Stent Mid Term	42
-------------------------	----

Ureteral Stent Long Term	43
--------------------------	----

PCN Catheters	48
---------------	----

Malecot	50
---------	----

Accessories

Bladder Evacuator	51
-------------------	----

Who we are

GBUK designs, develops and distributes millions of medical devices across multiple therapy areas including wound management, enteral feeding and patient handling. We are proud to be a top 20 supplier to the NHS, as well as supporting healthcare professionals and patients across the globe.

Our products are backed by comprehensive support, training and expert clinical teams. We are MDD and MDR compliant and hold ISO13485 Quality management system and ISO14001 Environmental management system.

Committed to environmental sustainability with self generated on-site electricity from our wind turbine and solar panels.

Our Mission:


GBUK, all subsidiaries and divisions are committed to improving patient outcomes and quality of life by delivering innovative, reliable, and easily accessible medical devices that empower healthcare professionals and serve global communities. As a trusted industry leader, we are committed to upholding the highest standards of regulatory compliance, ensuring safety and reliability in every product we offer, while creating lasting value for our customers and partners worldwide.



Access products

Guidewires

Guidewires provide reliable ureteral access and support smooth navigation throughout the urinary tract. They facilitate the placement, replacement, and exchange of instruments during a wide range of urological procedures, helping clinicians work with confidence and precision.




Key features

- Superior flexibility – promotes advancement through challenging or tortuous anatomy
- Enhanced durability – resists kinking and damage to deliver consistent, dependable performance
- Radiopaque – ensures clear visibility under fluoroscopy for accurate positioning

Access products

Striped Cobra Guidewire with Hydrophilic Tip


The hydrophilic tip enables smooth, atraumatic entry, while the striped shaft design delivers the pushability needed for controlled advancement. Constructed with nitinol and Teflon, the guidewire offers an optimal balance of flexibility, strength, and handling.



Features and benefits

- Hydrophilic tip for smooth atraumatic entry
- Constructed of Nitinol and Teflon
- Available in straight and angled tips

Diameter:	Length:	Tip:	Product code:
0.025"	260cm	Angled	MFCGW1425260
0.032"	150cm	Angled	MFCGW1432150
0.035"	150cm	Angled	MFCGW1435150
0.038"	150cm	Angled	MFCGW1438150
0.025"	260cm	Straight	MFCGW1325260
0.032"	150cm	Straight	MFCGW1332150
0.035"	150cm	Straight	MFCGW1335150
0.035"	260cm	Straight	MFCGW1335260
0.038"	150cm	Straight	MFCGW1338150



Straight




Angled

*Other sizes are available on request

Access products

Striped Guidewire

Medium stiffness provides stable pushability and controlled device passage, supporting confident advancement. The striped pattern enhances endoscopic visibility, enabling precise navigation during challenging cases. Constructed from Nitinol and Teflon, the guidewire delivers an ideal balance of strength, flexibility, and smooth handling.



Features and benefits

- The striped pattern enhances endoscopic visibility, enabling precise navigation during challenging cases
- Constructed of Nitinol and Teflon
- Available in straight and angled tips

Diameter:	Length:	Tip:	Product code:
0.035"	150cm	Angled	MFCGW1235150
0.035"	150cm	Straight	MFCGW1535150
0.038"	150cm	Straight	MFCGW1135150

*Other sizes are available on request



Straight




Angled

Access products

Nano Glide Hydrophilic Guidewire


The Nitinol and Teflon hydrophilic guidewire is engineered to streamline access by enabling rapid entry and smooth passage through obstructions, helping reduce overall procedure time.



Features and benefits

- Hydrophilic coated
- Constructed of Nitinol and Teflon
- Available in straight and angled tips

Diameter:	Length:	Tip:	Product code:
0.032"	150cm	Angled	MFCGW2432150
0.035"	150cm	Angled	MFCGW2435150
0.032"	150cm	Straight	MFCGW2332150
0.035"	150cm	Straight	MFCGW2335150
0.038"	150cm	Straight	MFCGW2338150



Straight




Angled

*Other sizes are available on request

Access products

PTFE Guidewire


Designed with a flexible tip and a moveable stainless-steel core, this guidewire offers control throughout the procedure. Its low-friction surface supports smooth passage, precise manoeuvrability, and reduced risk of kinking, even in challenging anatomy. Available in standard and super stiff options.



Features and benefits

- Moveable stainless-steel core
- Available in straight and J-tip types

Diameter:	Length:	Tip:	Stiffness:	Product code:
0.032"	150cm	J	Standard	MFCGW0232150
0.035"	150cm	J	Standard	MFCGW0235150
0.038"	150cm	J	Standard	MFCGW0238150
0.032"	150cm	Straight	Standard	MFCGW0132150
0.035"	150cm	Straight	Standard	MFCGW0135150
0.038"	150cm	Straight	Standard	MFCGW0138150
0.038"	80cm	J	Super stiff	MFCGW0338080
0.035"	150cm	Straight	Super stiff	MFCGW0435150
0.038"	150cm	Straight	Super stiff	MFCGW0438150
0.038"	260cm	Straight	Super stiff	MFCGW0438260



Straight




J-tip

*Other sizes are available on request

Access products

Ureteral Catheters

Used for drainage, retrograde pyelograms, atraumatic entry into the ureter, and for navigating tortuous ureters, this device supports smooth, controlled access during a range of urological procedures.




Key features

- Polyurethane with integrated Luer connection – provides smooth, secure, and reliable connectivity
- Atraumatic design – helps minimise patient discomfort when passing through the ureteral orifice
- Versatile functionality – suitable for retrograde pyelograms, drainage, and the administration of contrast medium
- Supplied with stylet

Access products

Ureteral Catheter Open End

Open-end design for guidewire insertion



Features and benefits

- Smooth advancement when navigating the ureter and accessing the renal pelvis
- Direct contrast injection capability – allows contrast to be delivered straight into the ureter or kidney during retrograde pyelography, enhancing diagnostic imaging
- Efficient drainage – the open end supports effective evacuation of urine or purulent material from the upper urinary tract


Diameter:	Length:	Coating:	Product code:
3Fr	70cm	N/A	MFUUC2203070
4Fr	70cm	N/A	MFUUC2204070
5Fr	70cm	N/A	MFUUC2205070
5Fr	70cm	Hydrophilic coating	MFUUC2205070H
6Fr	70cm	N/A	MFUUC2206070

*Other sizes are available on request

Access products

Ureteral Catheter Closed End

The closed-tip design reduces the risk of trauma to the ureteral mucosa during insertion, as there is no distal opening that could catch on tissue.



Features and benefits

- Prevents occlusion – side holes allow for effective drainage or infusion while minimising the risk of the tip becoming obstructed by debris or tissue.
- Uniform fluid distribution – Ideal for the even delivery of irrigating solutions or medications along the ureter


Diameter:	Length:	Coating:	Product code:
3Fr	70cm	N/A	MFUUC2103070
4Fr	70cm	N/A	MFUUC2104070

*Other tip types and sizes are available on request

Access products

Ureteral Catheter Cone Tip

The cone-shaped tip gently dilates the ureter as it advances, supporting easier passage through narrow or strictured areas.



Features and benefits

- Smooth navigation due to cone shaped tip
- Reduced trauma - Its tapered design helps minimise mucosal injury during insertion, promoting safer, more controlled access.
- Improved visualisation - The cone tip retains contrast media, enhancing fluoroscopic visibility for more precise assessment and navigation.


Diameter:	Length:	Coating:	Product code:
5Fr	70cm	N/A	MFUUC1305070

*Other tip types and sizes are available on request

Access products

Dual Lumen Catheter

The dual lumen design allows for precise guidewire insertion and efficient contrast media injection.



Features and benefits

- Smooth insertion – The smooth surface supports easy advancement while helping to minimise ureteral trauma.
- Effective ureteral dilation – A tapered tip design enables controlled and efficient dilation during access.
- Enhanced visualisation – Constructed from radiopaque material, the device provides clear fluoroscopic visibility for accurate placement.

Diameter:	Length:	Coating:	Product code:
10Fr	50cm	N/A	MFUDLC0110050
10Fr	50cm	Hydrophilic	MFUDLC0110050H

*Other sizes are available on request

Access products

Initial Puncture Needle

Used for percutaneous nephrostomy, cyst puncture, deep tissue access and all other percutaneous applications.



Features and benefits

- Compatible with all standard guidewire sizes
- Sharp tip for precise insertion and reduced tissue trauma
- Available in bevel, diamond and trocar tip
- Available in 2 or 3 parts
- Ergonomic hub design for improved control during insertion


Diameter:	Length:	Tip:	Part:	Product code:
18G	20cm	Bevel	2 part	MFCIPN1201820B

*Other sizes are available on request

Access products

Chiba Needle

A universal needle for percutaneous nephrostomy procedures, cytological biopsies, needle-aspirations, thoracentesis exploratory puncture and to inject contrast dyes.



Features and benefits

- Enables soft penetration and avoids the risk of vessel laceration.
- Ergonomic handling


Diameter:	Length:	Tip:	Product code:
16G	23cm	Bevel	MFCCN0101623
18G	20cm	Bevel	MFCCN0101820
22G	20cm	Bevel	MFCCN0102220

*Other sizes are available on request

Dilation

Urethral Dilator Set

Dilation of the urethra to establish and maintain access into the bladder for cystoscopy.



Features and benefits

- Smooth surface – facilitates easy introduction and minimising friction
- Atraumatic design – allowing for controlled and precise dilation
- Uniform taper – gradual dilation whilst reducing trauma to urethral tissues


Diameter:	Length:	Product code:
6-16Fr	40cm	MFUUD3161640
8-24Fr	40cm	MFUUD3182440

*Other sizes are available on request

Dilation

Urethral S Curve Dilator

Dilation of male urethra strictures. The curvature supports the natural anatomy of the male urethra, ensuring effective dilation.



Features and benefits

- Smooth surface – facilitates easy introduction and minimising friction
- Side hole for drainage
- Available individually or as a set

Dilator

Diameter:	Length:	Coating:	Product code:
22Fr	40cm	Hydrophilic	MFUSD0422040H
24Fr	40cm	Hydrophilic	MFUSD0424040H
26Fr	40cm	Hydrophilic	MFUSD0426040H

Dilator Set


Diameter:	Length:	Coating:	Product code:
8-20Fr	40cm	N/A	MFUSD0382040
8-20Fr	40cm	Hydrophilic	MFUSD0382040H

*Other sizes are available on request

Dilation

Ureteral Balloon Dilator

Used for transluminal dilation of ureteral strictures and for preparing the ureter prior to ureteroscopy or stone manipulation.



Features and benefits

- Radiopaque markers indicate the proximal and distal ends of the balloon for fluoroscopy placement
- Short tapered distal tip allows close approach to ureteral stone
- Made up of high strength material with excellent refolding characteristics
- Inflation pressure up to 16 ATM
- Provides effective and uniform dilation


Length:	Balloon diameter:	Balloon length:	Product code:
75cm	3.5mm	40mm	MFCBD01354075
75cm	4mm	40mm	MFCBD0144075
75cm	4mm	60mm	MFCBD0146075
75cm	5mm	40mm	MFCBD0154075
75cm	5mm	60mm	MFCBD0156075
75cm	6mm	40mm	MFCBD0164075
75cm	6mm	60mm	MFCBD0166075
75cm	7mm	60mm	MFCBD0176075

*Other sizes are available on request

Dilation

Ureteral Dilator Set

Dilation of the ureter prior to ureterorenoscopy, as well as step-by-step dilation of the intramural tunnel and ureter, or for precise stone manipulation.



Features and benefits

- Smooth surface – facilitates easy introduction and minimising friction
- Tapered tip – controlled insertion and reducing trauma during insertion
- Available with hydrophilic coating


Diameter:	Length:	Coating:	Product code:
6-16Fr	70cm	Hydrophilic	MFUUD3361670H
6-16Fr	70cm	N/A	MFUUD3061670
6-18Fr	70cm	N/A	MFUUD1261860

*Other sizes are available on request

Dilation

Ureteral Access Sheath - Braided

Used for ureteral dilation and to provide a working channel for the introduction of instruments during RIRS procedures.



Features and benefits

- Tapered tip for atraumatic insertion
- The braided tube has a Radiopaque marker band at the distal tip and dilator made up of radiopaque material for better (fluoroscopic) visualisation
- Available with hydrophilic coating

Size ID/OD:	Length:	Coating:	Product code:
9/11Fr	30cm	N/A	MFUUAS0409030
10/12Fr	28cm	Hydrophilic	MFUUAS0310028
10/12Fr	30cm	Hydrophilic	MFUUAS0310030
10/12Fr	35cm	Hydrophilic	MFUUAS0310035
10/12Fr	38cm	Hydrophilic	MFUUAS0310038
10/12Fr	45cm	Hydrophilic	MFUUAS0310045
10.7/12.7Fr	28cm	Hydrophilic	MFUUAS0310728
10.7/12.7Fr	35cm	Hydrophilic	MFUUAS0310735
11/13Fr	28cm	Hydrophilic	MFUUAS0311028
11/13Fr	45cm	Hydrophilic	MFUUAS0311045
12/14Fr	25cm	Hydrophilic	MFUUAS0312025
12/14Fr	28cm	Hydrophilic	MFUUAS0312028
12/14Fr	35cm	Hydrophilic	MFUUAS0312035
12/14Fr	38cm	Hydrophilic	MFUUAS0312038
12/14Fr	45cm	Hydrophilic	MFUUAS0312045
12/14Fr	38cm	Hydrophilic	MFUUAS0312538
14/16Fr	35cm	N/A	MFUUAS0414035
14/16Fr	45cm	N/A	MFUUAS0414045




Braided Hydrophilic
Ureteral Access Sheath

*Other sizes are available on request

Dilation

Ureteral Access Sheath - Non-Braided


Used for ureteral dilation and to provide a working channel for the introduction of instruments during RIRS procedures.



Features and benefits

- Tapered tip for atraumatic insertion
- The braided tube has a Radiopaque marker band at the distal tip and dilator made up of radiopaque material for better (fluoroscopic) visualisation
- Available with hydrophilic coating

Size ID/OD:	Length:	Coating:	Product code:
10/12Fr	35cm	Hydrophilic	MFUUAS0210035
11/13Fr	38cm	Hydrophilic	MFUUAS0211038
11/13Fr	45cm	Hydrophilic	MFUUAS0211045
11.5/13.5Fr	28cm	Hydrophilic	MFUUAS0211528
11.5/13.5Fr	35cm	Hydrophilic	MFUUAS0211535
12/14Fr	35cm	Hydrophilic	MFUUAS0212035
12/14Fr	45cm	Hydrophilic	MFUUAS0212045




Non-Braided Hydrophilic
Ureteral Access Sheath

*Other sizes are available on request

Dilation

Fascial Dilator

Used for dilation of the percutaneous tract, enabling safe insertion of surgical instruments.



Features and benefits

- Smooth Surface: Ensures easy introduction and reduces friction during dilation
- Uniform Taper: Allows for controlled and gradual dilation of fascial layers, minimising trauma
- Ergonomic Design: Ensures comfortable handling and precise control
- Available individually or as a set

Dilator

Diameter:	Length:	Coating:	Product code:
6Fr	20cm	Hydrophilic	MFCFD22006020H
8Fr	20cm	Hydrophilic	MFCFD22008020H
10Fr	20cm	Hydrophilic	MFCFD22010020H
12Fr	20cm	Hydrophilic	MFCFD22012020H
14Fr	20cm	Hydrophilic	MFCFD22014020H
16Fr	20cm	Hydrophilic	MFCFD22016020H
20Fr	20cm	Hydrophilic	MFCFD22020020H

Dilator Set

Diameter:	Length:	Coating:	Product code:
6-16Fr	20cm	Hydrophilic	MFCFD13061620H
8-16Fr	30cm	Hydrophilic	MFCFD13081630H

*Other sizes are available on request

Dilation

Amplatz Dilator

Used for dilation of the nephrostomy tract, enabling safe insertion and passage of surgical instruments.



Features and benefits

- Tapered tip design ensures smooth and atraumatic placement
- Anti-kink technology: Firm, kink-resistant material ensures reliable dilation without bending or collapsing
- Available individually, with sheath or as a set

Amplatz Dilator

Diameter:	Length:	Product code:
12Fr	35cm	MFUAD011235
18Fr	35cm	MFUAD011835

Amplatz Sheath

Diameter:	Length:	Product code:
24/28Fr	17cm	MFUAS01282417
26/30Fr	17cm	MFUAS01302617
28/32Fr	17cm	MFUAS01322817

Amplatz Dilator with Sheath

Diameter:	Length:	Product code:
20Fr	30cm	MFUASD012030
30Fr	30cm	MFUASD013030

Amplatz Renal Dilator Set


Diameter:	Length:	Product code:
6-30Fr	30cm	MFUADS0206030LLT
6-30Fr	35cm	MFUADS0206035LLT

*Other sizes are available on request

Dilation

Nephrostomy Balloon Dilator

Used to dilate the musculofascial tract, and renal capsule establishing a secure Nephrostomy tract during Percutaneous Nephrolithotomy (PCNL).



Features and benefits

- Anti-kink design minimizes kinking and bending, reducing procedural complications
- Tapered tip ensures smooth tract dilation while minimising trauma to surrounding tissues
- Radiopaque markers enhance visualisation of balloon position under fluoroscopy


Length:	Balloon diameter:	Balloon length:	Inflation device	Product code:
55cm	10mm	12cm	No	MFCNB01001255
55cm	10mm	12cm	Yes	MFCNB01001255ID

*Other sizes are available on request

Stone management

G-Paw


Combining the benefits of a stone grasper and a stone basket. It allows for easy grasp, release and removal of stones.



Features and benefits

- Nitinol wire, known for its shape retention, enhances manoeuvrability and precision during procedure
- Tipless design permits atraumatic capture of stone

Balloon diameter:	Basket diameter:	Length:	Number of wires:	Product code:
1.9Fr	9mm	120cm	3	MFUGP91193120F
1.9Fr	12mm	120cm	3	MFUGP92193120F
2.2Fr	9mm	120cm	3	MFUGP91223120F
2.2Fr	12mm	120cm	3	MFUGP92223120F




*Other sizes are available on request

Stone management

Tipless Stone Basket

Used for retrieval of kidney and ureteric stones.




Key features

- Nitinol wire, known for its shape retention, enhances manoeuvrability and precision during procedures
- Handle permits complete disassembly
- Tipless design permits atraumatic capture of stone
- Available in multiple handle types

Stone management

Tipless Stone Basket - Balloon

A flexible basket that expands into a balloon shape, surrounding the stones for easy extraction and minimising trauma to surrounding tissues.



Features and benefits


- Nitinol wire, known for its shape retention, enhances manoeuvrability and precision during procedure
- 4 wire basket
- Available in variable basket diameters

Diameter:	Basket diameter:	Length:	Number of wires:	Handle:	Product code:
1.7Fr	12mm	120cm	4	Fish	MFUSB18174120F
1.7Fr	15mm	120cm	4	Fish	MFUSB28174120F
1.9Fr	12mm	90cm	4	Fish	MFUSB18194090F
1.9Fr	12mm	120cm	4	Fish	MFUSB18194120F
1.9Fr	15mm	120cm	4	Fish	MFUSB28194120F
2.2Fr	12mm	90cm	4	Fish	MFUSB18224090F
2.2Fr	12mm	120cm	4	Fish	MFUSB18224120F
2.2Fr	15mm	120cm	4	Fish	MFUSB28224120F
2.2Fr	12mm	120cm	4	Scissor	MFUSB18224120S
3.0Fr	12mm	120cm	4	Fish	MFUSB18304120F
3.0Fr	15mm	120cm	4	Fish	MFUSB28304120F
4.0Fr	15mm	70cm	4	Fish	MFUSB28404070F

Stone management

Tipless Stone Basket - Segura

Helical or spiral geometry that wraps around stones to ensure a firm grasp for retrieval in difficult to reach areas.



Features and benefits


- Nitinol wire, known for its shape retention, enhances manoeuvrability and precision during procedure
- 4 wire basket
- Available in variable basket diameters

Diameter:	Basket diameter:	Length:	Number of wires:	Handle:	Product code:
1.9Fr	15mm	90cm	4	Fish	MFUSB24194090F
1.9Fr	15mm	120cm	4	Fish	MFUSB24194120F
1.9Fr	12mm	120cm	4	Scissor	MFUSB14194120S
2.4Fr	12mm	120cm	4	Scissor	MFUSB14244120S
3.0Fr	15mm	90cm	4	Scissor	MFUSB54304090S
3.0Fr	15mm	120cm	4	Scissor	MFUSB54304120S
4.0Fr	9mm	90cm	4	Fish	MFUSB04404090F
4.0Fr	9mm	120cm	4	Fish	MFUSB04404120F

Stone management

Tipless Stone Basket - Helical

A multi-wire basket that adapts to the shape of the stone, offering excellent flexibility and enhanced manoeuvrability for efficient retrieval.



Features and benefits


- Nitinol wire, known for its shape retention, enhances manoeuvrability and precision during procedure
- 4 wire basket
- Available in variable basket diameters

Diameter:	Basket diameter:	Length:	Number of wires:	Handle:	Product code:
1.7Fr	12mm	120cm	4	Fish	MFUSB16174120F
1.9Fr	12mm	90cm	4	Scissor	MFUSB46194090S
1.9Fr	12mm	120cm	4	Scissor	MFUSB46194120S
2.2Fr	12mm	90cm	4	Scissor	MFUSB46224090S
2.2Fr	12mm	120cm	4	Scissor	MFUSB46224120S
3.0Fr	15mm	120cm	4	Fish	MFUSB26304120F

Stone management

Tipless Stone Basket - X-Circle

Provides a unique net to capture and extract intact stones and prevent migration.



Features and benefits


- Nitinol wire, known for its shape retention, enhances manoeuvrability and precision during procedure
- 12, 16 wire basket
- Available in variable basket diameters

Diameter:	Basket diameter:	Length:	Number of wires:	Handle:	Product code:
1.9Fr	12mm	90cm	16	Scissor	MFUSX241916090S
1.9Fr	12mm	120cm	16	Fish	MFUSX241916120F
1.9Fr	15mm	120cm	16	Fish	MFUSX251916120F
1.9Fr	18mm	120cm	16	Fish	MFUSX261916120F
1.9Fr	12mm	120cm	16	Fish	MFUSX241916120S
2.2Fr	12mm	120cm	12	Scissor	MFUSX312212120S
2.2Fr	12mm	90cm	16	Scissor	MFUSX242216090S
2.2Fr	18mm	120cm	16	Fish	MFUSX262216120F
2.2Fr	12mm	120cm	16	Scissor	MFUSX242216120S
2.5Fr	15mm	120cm	16	Fish	MFUSX252516120F
3.0Fr	12mm	90cm	12	Scissor	MFUSX213012090S
3.0Fr	12mm	120cm	16	Scissor	MFUSX243016120S

Stone management

Tipped Stone Basket

Features a tapered distal tip that aids in navigation and controlled positioning for efficient retrieval.




Key features

- Nitinol wire, known for its shape retention, enhances manoeuvrability and precision during procedures
- 4 wire configuration
- Handle permits complete disassembly
- Available in multiple handle types

Stone management

Tipped Stone Basket - Helical

A multi-wire basket that adapts to the shape of the stone, offering excellent flexibility and enhanced manoeuvrability for efficient retrieval.



Features and benefits


- Nitinol wire, known for its shape retention, enhances manoeuvrability and precision during procedure
- 4 wire basket
- Available in variable basket diameters

Diameter:	Basket diameter:	Length:	Number of wires:	Handle:	Product code:
2.2Fr	9mm	120cm	4	Scissor	MFUSB05224120S
2.5Fr	9mm	90cm	4	Fish	MFUSB05254090F
2.5Fr	9mm	120cm	4	Fish	MFUSB05254120F
2.5Fr	12mm	120cm	4	Fish	MFUSB15254120F
2.5Fr	12mm	120cm	4	Fish	MFUSB45254120F
2.5Fr	9mm	120cm	4	Scissor	MFUSB01254120S
3.0Fr	9mm	90cm	4	Fish	MFUSB05304090F
3.0Fr	9mm	120cm	4	Fish	MFUSB05304120F
3.0Fr	12mm	120cm	4	Fish	MFUSB15304120F
3.0Fr	15mm	120cm	4	Fish	MFUSB25304120F
3.0Fr	12mm	120cm	4	Fish	MFUSB45304120F
4.0Fr	12mm	90cm	4	Scissor	MFUSB11404090S
4.0Fr	12mm	120cm	4	Scissor	MFUSB11404120S

Stone management

Tipped Stone Basket - Segura

Designed for reliable stone capture, the Segura tip stone basket features a tapered distal tip that aids in navigation and controlled positioning for efficient retrieval.



Features and benefits


- Nitinol wire, known for its shape retention, enhances manoeuvrability and precision during procedure
- 4 wire basket
- Available in variable basket diameters

Diameter:	Basket diameter:	Length:	Number of wires:	Handle:	Product code:
2.5Fr	15mm	90cm	4	Fish	MFUSB23254090F
2.5Fr	15mm	120cm	4	Fish	MFUSB23254120F
3.0Fr	15mm	90cm	4	Fish	MFUSB23304090F
3.0Fr	12mm	120cm	4	Fish	MFUSB13304120F
3.0Fr	15mm	120cm	4	Fish	MFUSB23304120F

Stone management

Perk Basket

Used for kidney stone removal.



Features and benefits

- Available as G-paw, Segura or X-circle tip

Tipless Segura basket option

Shaft

Handle assembly

Diameter:	Basket Diameter:	Length:	Tip:	Product code:
8Fr	17mm	38cm	Segura	MFUPB01083817
10Fr	17mm	38cm	Segura	MFUPB01103817
10Fr	N/A	38cm	G-Paw	MFUPB10010038
10Fr	N/A	38cm	X-Circle	MFUPB11010038
12Fr	N/A	38cm	X-Circle	MFUPB11012038

*Other sizes are available on request

Basket options



G-Paw



Segura




X-Circle

Stone management

Stone Grasper

Used in the facilitation of endoscopic manipulation or removal of stones and retrieval of foreign object.



Features and benefits

- Tri-prong, atraumatic design, minimises tissue damage during procedures
- Maximum grasping efficiency, ensures secure handling and removal of stones and foreign objects


Diameter:	Basket Diameter:	Length:	Handle:	Product code:
2.5Fr	15mm	120cm	Fish	MFUSG0325120F
3.2Fr	12mm	115cm	Scissor	MFUSG1232115S

*Other sizes are available on request

Stone management

Stone Stop

Prevents stone migration to the kidney during lithotripsy.




Diameter:	Length:	Handle:	Product code:
3Fr	120cm	Clip	MFUSS0110120C
3Fr	120cm	Push	MFUSS0110120P

*Other sizes are available on request

Stone management

URS Forceps

Used for grasping, manipulating and extracting calculi from ureteral and renal pelvis.



Features and benefits

- Different handles in URS forceps enables choice to the surgeons in preferring different ergonomics that suits their style of working
- The rigid handle is the standard model, while the flexible handle variant is designed to access multiple calyces, helping minimise the need for multiple punctures during complete stone clearance


Diameter:	Length:	Tip:	Product code:
3Fr	65cm	Aligator	MFUF043065

*Other sizes and tips are available on request

Drainage

Ureteral Double J Stent

Used for internal drainage from the ureteropelvic junction to the bladder.




Key features

- Highly radiopaque for better visualisation and has an extremely low encrustation tendency

Drainage

Ureteral Double J Stent - Medium Term

For mid-term use (30–90 days) providing steady drainage support with enhanced flexibility, excellent visibility, and a smooth patient experience.



Features and benefits

- Made of Carbothane
- Radiopaque for easy visualisation
- Available individually or as a set
- Available with hydrophilic coating

Medium Term Stent

Diameter:	Length:	Tip configuration:	Multiloop:	Product code:
3.5Fr	8cm	Both ends open	Yes	MFDJA7603508
3.5Fr	22cm	Both ends open	Yes	MFDJA7603522
4.7Fr	16cm	Both ends open	Yes	MFDJA7604716

Medium Term Stent Set


Diameter:	Length:	Tip configuration:	Tether:	Clamp:	Guidewire:	Product code:
3.0Fr	8cm	Both ends open	Yes	Yes	Yes	MFDJD5403008
3.5Fr	8cm	Both ends open	Yes	Yes	Yes	MFDJD5403508
3.5Fr	12cm	Both ends open	Yes	Yes	Yes	MFDJD5403512
3.5Fr	14cm	Both ends open	Yes	Yes	Yes	MFDJD5403514
3.5Fr	16cm	Both ends open	Yes	Yes	Yes	MFDJD5403516
6.0Fr	12cm	Both ends open	Yes	Yes	Yes	MFDJD5406012

*Other sizes are available on request

Drainage

Ureteral Double J Stent - Long Term

Engineered for long-term ureteral drainage (up to 360 days). They are biocompatible, flexible, and radiopaque, providing lasting performance with minimal encrustation risk. Silicone is considered better suited for reduced encrustation in long-term DJ (Double-J) stents due to a combination of material properties that limit mineral deposition, bacterial adherence, and biofilm formation.



Features and benefits

- Made of silicone
- Radiopaque for easy visualisation
- Available individually or as a set
- Available with hydrophilic coating

Drainage

Ureteral Double J Stent - Long Term

Long Term Stent

Diameter:	Length:	Tip configuration:	Coating:	Product code:
4.5Fr	22cm	Both ends open	Hydrophilic	MFDJA2404522H
5.0Fr	24cm	Both ends open		MFDJA2405024
6.0Fr	22cm	Both ends open		MFDJA2406022
6.0Fr	22cm	Both ends open	Hydrophilic	MFDJA2406022H
6.0Fr	24cm	Both ends open		MFDJA2406024
6.0Fr	24cm	Both ends open	Hydrophilic	MFDJA2406024H
6.0Fr	26cm	Both ends open		MFDJA2406026
6.0Fr	26cm	Both ends open	Hydrophilic	MFDJA2406026H
6.0Fr	28cm	Both ends open		MFDJA2406028
6.0Fr	28cm	Both ends open	Hydrophilic	MFDJA2406028H
7.0Fr	22cm	Both ends open		MFDJA2407022
7.0Fr	24cm	Both ends open	Hydrophilic	MFDJA2407024H
8.0Fr	22cm	Both ends open	Hydrophilic	MFDJA2408022H
8.0Fr	24cm	Both ends open	Hydrophilic	MFDJA2408024H
8.0Fr	26cm	Both ends open		MFDJA2408026
8.0Fr	26cm	Both ends open	Hydrophilic	MFDJA2408026H
8.0Fr	28cm	Both ends open		MFDJA2408028

Diameter:	Length:	Tip configuration:	Coating:	Multiloop:	Product code:
4.7Fr	22cm	Both ends open		Yes	MFDJA2604722
4.8Fr	22cm	Both ends open	Hydrophilic	Yes	MFDJA2604822H
5.0Fr	22cm	Both ends open		Yes	MFDJA2605022
6.0Fr	22cm	Both ends open		Yes	MFDJA2606022
6.0Fr	22cm	Both ends open	Hydrophilic	Yes	MFDJA2606022H
7.0Fr	22cm	Both ends open		Yes	MFDJA2607022
7.0Fr	22cm	Both ends open	Hydrophilic	Yes	MFDJA2607022H

Drainage

Ureteral Double J Stent - Long Term

Long Term Stent with Suture

Diameter:	Length:	Tip configuration:	Coating:	Product code:
4.5Fr	24cm	Both ends open	Hydrophilic	MFDJC2404524H
4.5Fr	26cm	Both ends open	Hydrophilic	MFDJC2404526H
4.8Fr	22cm	Both ends open	Hydrophilic	MFDJC2404822H
4.8Fr	24cm	Both ends open	Hydrophilic	MFDJC2404824H
4.8Fr	26cm	Both ends open	Hydrophilic	MFDJC2404826H
5.0Fr	26cm	Both ends open	Hydrophilic	MFDJC2405026H
6.0Fr	20cm	Both ends open	Hydrophilic	MFDJC2406020H
6.0Fr	22cm	Both ends open	Hydrophilic	MFDJC2406022H
6.0Fr	24cm	Both ends open	Hydrophilic	MFDJC2406024H
6.0Fr	26cm	Both ends open		MFDJC2406026
6.0Fr	26cm	Both ends open	Hydrophilic	MFDJC2406026H
6.0Fr	28cm	Both ends open	Hydrophilic	MFDJC2406028H
6.0Fr	30cm	Both ends open	Hydrophilic	MFDJC2406030H
7.0Fr	22cm	Both ends open	Hydrophilic	MFDJC2407022H
7.0Fr	24cm	Both ends open	Hydrophilic	MFDJC2407024H
7.0Fr	26cm	Both ends open	Hydrophilic	MFDJC2407026H
8.0Fr	22cm	Both ends open	Hydrophilic	MFDJC2408022H
8.0Fr	24cm	Both ends open	Hydrophilic	MFDJC2408024H
8.0Fr	26cm	Both ends open	Hydrophilic	MFDJC2408026H
6.0Fr	12cm	Closed end	Hydrophilic	MFDJC2106012H
6.0Fr	24cm	Closed end	Hydrophilic	MFDJC2106024H
6.0Fr	26cm	Closed end	Hydrophilic	MFDJC2106026H
6.0Fr	28cm	Closed end	Hydrophilic	MFDJC2106028H

Diameter:	Length:	Tip configuration:	Coating:	Multiloop:	Product code:
6.0Fr	22cm	Both ends open	Hydrophilic	Yes	MFDJC2606022H
7.0Fr	22cm	Both ends open	Hydrophilic	Yes	MFDJC2607022H
8.0Fr	22cm	Both ends open	Hydrophilic	Yes	MFDJC2608022H
4.5Fr	22cm	Closed end	Hydrophilic	Yes	MFDJC2304522H

Drainage

Ureteral Double J Stent - Long Term

Long Term Stent including Suture and Guidewire

Diameter:	Length:	Tip configuration:	Coating:	Product code:
4.5Fr	12cm	Both ends open		MFDJD2404512
4.5Fr	14cm	Both ends open		MFDJD2404514
4.7Fr	16cm	Both ends open		MFDJD2404716
4.7Fr	18cm	Both ends open		MFDJD2404718
4.7Fr	22cm	Both ends open		MFDJD2404722
4.7Fr	26cm	Both ends open		MFDJD2404726
4.8Fr	24cm	Both ends open	Hydrophilic	MFDJD2404824H
4.8Fr	26cm	Both ends open	Hydrophilic	MFDJD2404826H
4.8Fr	28cm	Both ends open		MFDJD2404828
4.8Fr	28cm	Both ends open	Hydrophilic	MFDJD2404828H
5.0Fr	26cm	Both ends open		MFDJD2405026
6.0Fr	12cm	Both ends open		MFDJD2406012
6.0Fr	14cm	Both ends open		MFDJD2406014
6.0Fr	16cm	Both ends open		MFDJD2406016
6.0Fr	22cm	Both ends open	Hydrophilic	MFDJD2406022H
6.0Fr	24cm	Both ends open	Hydrophilic	MFDJD2406024H
6.0Fr	26cm	Both ends open		MFDJD2406026
6.0Fr	26cm	Both ends open	Hydrophilic	MFDJD2406026H
7.0Fr	24cm	Both ends open	Hydrophilic	MFDJD2407024H
7.0Fr	26cm	Both ends open	Hydrophilic	MFDJD2407026H
7.0Fr	28cm	Both ends open	Hydrophilic	MFDJD2407028H
8.0Fr	24cm	Both ends open	Hydrophilic	MFDJD2408024H
8.0Fr	26cm	Both ends open	Hydrophilic	MFDJD2408026H
8.0Fr	28cm	Both ends open	Hydrophilic	MFDJD2408028H

Drainage

Ureteral Double J Stent - Long Term

Long Term Stent including Suture and Guidewire


Diameter:	Length:	Tip configuration:	Coating:	Multiloop:	Product code:
4.8Fr	22cm	Both ends open	Hydrophilic	Yes	MFDJD2604822H
5.0Fr	22cm	Both ends open		Yes	MFDJD2605022
6.0Fr	22cm	Both ends open		Yes	MFDJD2606022
6.0Fr	22cm	Both ends open	Hydrophilic	Yes	MFDJD2606022H
7.0Fr	22cm	Both ends open		Yes	MFDJD2607022
8.0Fr	22cm	Both ends open		Yes	MFDJD2608022
8.0Fr	22cm	Closed End		Yes	MFDJD2308022

*Other sizes are available on request

Drainage

PCN Catheter

Used for temporary or permanent drainage and irrigation of urine from the kidney by direct puncture, in the presence of distal obstruction.



Features and benefits

- Soft material variant for patient comfort
- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- Hydrophilic coating allows easy insertion
- Drainage efficiency

PCN Catheter Set contains:

- PCN Catheter
- Guidewire
- IP Needle
- Fascial dilators
- Scalpel blade

PCN Catheter

Diameter:	Length:	Coating:	Product code:
6Fr	30cm	Hydrophilic	MFCPC1606030H

PCN Catheter with Suture

Diameter:	Length:	Coating:	Product code:
6Fr	30cm	Hydrophilic	MFCPC2206030H
7Fr	30cm	Hydrophilic	MFCPC2207030H
8Fr	30cm	Hydrophilic	MFCPC2208030H
9Fr	30cm	Hydrophilic	MFCPC2209030H
10Fr	35cm	Hydrophilic	MFCPC2210035H
12Fr	35cm	Hydrophilic	MFCPC2212035H
14Fr	35cm	Hydrophilic	MFCPC2214035H

Drainage

PCN Catheter

PCN Catheter Set

Diameter:	Length:	Coating:	Product code:
10Fr	30cm	Hydrophilic	MFCPC4310030H
12Fr	35cm	Hydrophilic	MFCPC4312035H
14Fr	35cm	Hydrophilic	MFCPC4314035H

PCN/Abscess Drainage Catheter

Diameter:	Length:	Product code:
8Fr	25cm	MFRPC1308025
10Fr	25cm	MFRPC1310025
12Fr	25cm	MFRPC1312025

PCN Catheter with Locking Pigtail

Diameter:	Length:	Coating:	Tip configuration:	Product code:
6Fr	40cm	Hydrophilic	Locking pigtail	MFCPC3706400H
8fr	40cm	Hydrophilic	Locking pigtail	MFCPC3708040H
10Fr	40cm	Hydrophilic	Locking pigtail	MFCPC3701040H

PCN Catheter with Needle

Diameter:	Length:	Coating:	Product code:
9Fr	30cm	Hydrophilic	MFCPC1309030H
10Fr	35cm	Hydrophilic	MFCPC1310035H

PCN Catheter Set with Suture


Diameter:	Length:	Coating:	Product code:
10Fr	35cm	Hydrophilic	MFCPC6210035H
12Fr	35cm	Hydrophilic	MFCPC6212035H
14Fr	35cm	Hydrophilic	MFCPC6214035H
16Fr	35cm	Hydrophilic	MFCPC6216035H

*Other sizes are available on request

Drainage

Malecot Catheter

The Malecot flower design is ideally suited for the drainage of thick viscous fluids. Used for percutaneous placement of a malecot catheter in the renal pelvis for nephrostomy drainage.



Features and benefits

- Integral flower design matches with the catheter size when closed for less trauma during insertion
- Available in a Set

Malecot Catheter Set contains:

- Malecot Catheter with Pusher
- Guidewire
- IP Neele
- Fascial dilators
- Scalpel blade
- Urine bag connector
- 3 way stopcock
- Small clamp

Malecot Catheter

Diameter:	Length:	Product code:
10Fr	35cm	MFCMC2110035
12Fr	35cm	MFCMC2112035
14Fr	35cm	MFCMC2114035
16Fr	35cm	MFCMC2116035

Malecot Catheter Set

Diameter:	Length:	Product code:
10Fr	35cm	MFCMC2510035
12Fr	35cm	MFCMC2512035
14Fr	35cm	MFCMC2514035
16Fr	35cm	MFCMC2516035
18Fr	30cm	MFCMC6518030

*Other sizes are available on request

Accessories

Bladder Evacuator

Used to perform dual irrigating and tissue straining/collecting functions during transurethral bladder and prostate surgery. Compatible with Storz, Wolf and Olympus.



Features and benefits

- Self-air-purging design allows efficient and rapid tissue evacuation
- Adaptable connection end provides an easy, watertight connection to scope
- Airtight seal prevents risk of fluid or tissue loss and maintains suction pressure

Description:

Product code:

Bladder Evacuator

MFBE001

GBUK Group Ltd.

Woodland House, Blackwood Hall Business Park,
North Duffield, Selby, North Yorkshire,
YO8 5DD, UK

T +44 (0) 1757 288 587

E info@gbukgroup.com

W www.gbukgroup.com



GBUK
ENHANCING PATIENT CARE