

Compatibility Statement

Sling products within the Care & Independence GLOVE SLING™ and COMFORT VALUE™ ranges are deemed compatible with the hoist models from the manufacturers' listed below. In the event of any doubt or if a specific hoist does not appear in the list, please contact us to confirm compatibility *before* using the Care & Independence sling.

NOTE: GLOVE SLING™ products are manufactured with a loop strap configuration as standard but a clip variant may be available upon request. Please contact us for further details.

MANUFACTURER	Track Hoists	Mobile Hoists
Arjo	Maxi Sky-2, 2-plus, 1000 & 440	Maxi -500, Move, Twin, Minstrel & Tenor
Aldersly		Bariatric Hoist
BHM	Voyager, V4, Maxi	
Chiltern	Wispa 100, 200, 300, Wispa Lite	Wispa
Ergolet	Luna	
Gulldman	Dan-Hoist/BASIC. GH2, GH3.	GL5
Handicare	Rise Atlas, basic track hoists	Carina, Eva, Vega
Handi-Move		Victor
Helping Hand	Skyframe	Swan
Homecraft		Mobile hoist
Horcher	Diana* & Lexa <i>*Excludes GHUG Mk1</i>	
Huntleigh		Porta 127, 150, 155, 175 and 225
Invacare	Robin, 800	Birdie, Alto, Compact, Reliant and Saxo
Joerns/Oxford	Voyager 420, 550 & 800 Voyager 420, 550 & 800 Elara, Voyager	Advance, Presence, Mini, Midi & Maxi
Littech	Gemini	
Liko	Likorail, Multirail & Ultra Twin	Golvo, Viking and Uno
Macworth	Essential. M Series	
Molift	Air, Nomad	Smart, Mover & Partner
Monarch	Track hoist	
Morecare	Swift Lifter	
Nopac	Windsor, Royal	
Prism Medical	Freeway Transactive, Prism CP200	Freeway M Series, S Series, A Series
Smith & Nephew		Liftmaster
Southern Care	H130, H132, H190, HS200, EasyHoist	
Symetrik	Skyframe	
Unihoist Ltd		University Scholar, Graduate, Tutor and Master
Test Valley		Nova
Tri-Care	Freeway mobile, track and portable hoist	
Waverly Glen	Transactive, Transportable track hoist	
Wessex	Travelmaster	