



Enhancing patient care

Nutricare[®] Gastrostomy Product Range

Our Mission


GBUK, all subsidiaries and divisions are committed to improving patient outcomes and quality of life by delivering innovative, reliable, and easily accessible medical devices that empower healthcare professionals and serve global communities. As a trusted industry leader, we are committed to upholding the highest standards of regulatory compliance, ensuring safety and reliability in every product we offer, while creating lasting value for our customers and partners worldwide.

Contents

Nutricare® compact balloon gastrostomy tube	4
Nutricare® low profile balloon gastrostomy button kit	6
Nutricare® feed sets for Nutricare low profile balloon gastrostomy button	12
Nutricare® stoma measuring device	13
GBUK Home Hub	14

Nutricare® compact balloon gastrostomy tube

Introducing a new generation of gastrostomy feeding tube - designed for patients who fall between traditional options. When a low-profile button isn't appropriate and a longer gastrostomy tube is too cumbersome, our hybrid design offers the perfect balance of comfort, practicality, and discretion.



Nutricare® compact balloon gastrostomy tube

The Nutricare compact gastrostomy tube is intended to provide feed, medication, and/or hydration directly into the stomach, through the abdominal wall. It may also be used for gastric decompression. It is intended for placement into an established gastrostomy stoma tract, as well as for initial placement into a newly established gastrostomy stoma tract with the use of gastropexy sutures.



Features and benefits

- **Durable balloon**, so water loss is reduced due to gas barrier design, allowing for monthly water changes
- **Tapered tip** for ease of insertion
- **Radiopaque line** for visibility on X-ray
- **Variable fill balloon** to tailor balloon fill for better patient fit
- **Apple shaped balloon** for an effective seal at the stoma tract
- **Made of soft silicone** for patient comfort
- **Right angle external retention device** allows the tube to be folded securely against the skin creating a low profile and discreet device
- **Thoughtful external retention device design** allows for air circulation

Size	Description	UOI	Code	NHSSC
14Fr	Nutricare compact balloon gastrostomy tube	1	GBCGT14	FAR85080
16Fr	Nutricare compact balloon gastrostomy tube	1	GBCGT16	FAR85081
18Fr	Nutricare compact balloon gastrostomy tube	1	GBCGT18	FAR85082
20Fr	Nutricare compact balloon gastrostomy tube	1	GBCGT20	FAR85084
24Fr	Nutricare compact balloon gastrostomy tube	1	GBCGT24	FAR85083

Nutricare® low profile balloon gastrostomy button kit

The Nutricare Button is intended to provide feed, medication, and/or hydration directly into the stomach, through the abdominal wall. It may also be used for gastric decompression. It is intended for placement into an established gastrostomy stoma tract, as well as for initial placement into a newly established Gastrostomy stoma tract with the use of gastropexy sutures.



Tapered tip



Apple shaped balloon



Soft silicone material

Features and benefits


- **Extra small footprint** to allow for air circulation and comfort
- **Apple shaped balloon** for an effective seal at the stoma tract
- **Durable balloon**, so water loss is reduced due to gas barrier design, allowing for monthly water changes
- **Size specific introducer** to aid insertion if required
The tip of the introducer corresponds with the button shaft length, and will not protrude past the tube tip
- **Tapered tip** for ease of insertion
- **Variable fill balloon** to tailor balloon fill for better patient fit
- **Soft silicone** for patient comfort
- **Radiopaque line** for visibility on X-ray
- **Convenient** kit includes three feed sets, a size specific introducer, gauze and balloon inflation syringe
- **Clearly marked** balloon inflation port clearly marked with recommended fill volume, for convenience



Size specific introducer

Nutricare® low profile balloon gastrostomy button kit


The Nutricare Button is one of the smallest and most discreet gastrostomy buttons compared to major market alternatives.



Button	Skin contact area		Overall profile area	
Market alternative 1	291mm ²		438mm ²	
Market alternative 2	213mm ²		377mm ²	
Nutricare	193mm ²		299mm ²	

Nutricare® low profile balloon gastrostomy button kit

Each kit contains:




1 x 10ml luer slip balloon syringe




5 x gauze




1 x size specific introducer
(per button size)




1 x Nutricare™ low profile
balloon gastrostomy button



1 x 90mm
straight feed set



1 x 280mm
right angle feed set



1 x 280mm
straight feed set

Nutricare® low profile balloon gastrostomy button kit codes

Size	Description	UOI	Code	NHSSC
12Fr x 0.8cm	Low profile balloon gastrostomy button kit	1	GBX1208	FAN85372
12Fr x 1.0cm	Low profile balloon gastrostomy button kit	1	GBX1210	FAN85399
12Fr x 1.2cm	Low profile balloon gastrostomy button kit	1	GBX1212	FAN85355
12Fr x 1.5cm	Low profile balloon gastrostomy button kit	1	GBX1215	FAN85354
12Fr x 1.7cm	Low profile balloon gastrostomy button kit	1	GBX1217	FAN85398
12Fr x 2.0cm	Low profile balloon gastrostomy button kit	1	GBX1220	FAN85397
12Fr x 2.5cm	Low profile balloon gastrostomy button kit	1	GBX1225	FAN85396
12Fr x 3.0cm	Low profile balloon gastrostomy button kit	1	GBX1230	FAN85353
14Fr x 0.8cm	Low profile balloon gastrostomy button kit	1	GBX1408	FAN85352
14Fr x 1.0cm	Low profile balloon gastrostomy button kit	1	GBX1410	FAN85391
14Fr x 1.2cm	Low profile balloon gastrostomy button kit	1	GBX1412	FAN85390
14Fr x 1.5cm	Low profile balloon gastrostomy button kit	1	GBX1415	FAN85389
14Fr x 1.7cm	Low profile balloon gastrostomy button kit	1	GBX1417	FAN85388
14Fr x 2.0cm	Low profile balloon gastrostomy button kit	1	GBX1420	FAN85349
14Fr x 2.3cm	Low profile balloon gastrostomy button kit	1	GBX1423	FAN85385
14Fr x 2.5cm	Low profile balloon gastrostomy button kit	1	GBX1425	FAN85367
14Fr x 2.7cm	Low profile balloon gastrostomy button kit	1	GBX1427	FAN85408
14Fr x 3.0cm	Low profile balloon gastrostomy button kit	1	GBX1430	FAN85384
14Fr x 3.5cm	Low profile balloon gastrostomy button kit	1	GBX1435	FAN85366
14Fr x 4.0cm	Low profile balloon gastrostomy button kit	1	GBX1440	FAN85409
14Fr x 4.5cm	Low profile balloon gastrostomy button kit	1	GBX1445	FAN85369
16Fr x 0.8cm	Low profile balloon gastrostomy button kit	1	GBX1608	FAN85392
16Fr x 1.0cm	Low profile balloon gastrostomy button kit	1	GBX1610	FAN85370
16Fr x 1.2cm	Low profile balloon gastrostomy button kit	1	GBX1612	FAN85410
16Fr x 1.5cm	Low profile balloon gastrostomy button kit	1	GBX1615	FAN85411

Nutricare® low profile balloon gastrostomy button kit codes

Size	Description	UOI	Code	NHSSC
16Fr x 1.7cm	Low profile balloon gastrostomy button kit	1	GBX1617	FAN85350
16Fr x 2.0cm	Low profile balloon gastrostomy button kit	1	GBX1620	FAN85349
16Fr x 2.3cm	Low profile balloon gastrostomy button kit	1	GBX1623	FAN85385
16Fr x 2.5cm	Low profile balloon gastrostomy button kit	1	GBX1625	FAN85367
16Fr x 2.7cm	Low profile balloon gastrostomy button kit	1	GBX1627	FAN85357
16Fr x 3.0cm	Low profile balloon gastrostomy button kit	1	GBX1630	FAN85341
16Fr x 3.5cm	Low profile balloon gastrostomy button kit	1	GBX1635	FAN85400
16Fr x 4.0cm	Low profile balloon gastrostomy button kit	1	GBX1640	FAN85358
16Fr x 4.5cm	Low profile balloon gastrostomy button kit	1	GBX1645	FAN85375
18Fr x 0.8cm	Low profile balloon gastrostomy button kit	1	GBX1808	FAN85359
18Fr x 1.0cm	Low profile balloon gastrostomy button kit	1	GBX1810	FAN85360
18Fr x 1.2cm	Low profile balloon gastrostomy button kit	1	GBX1812	FAN85361
18Fr x 1.5cm	Low profile balloon gastrostomy button kit	1	GBX1815	FAN85342
18Fr x 1.7cm	Low profile balloon gastrostomy button kit	1	GBX1817	FAN85343
18Fr x 2.0cm	Low profile balloon gastrostomy button kit	1	GBX1820	FAN85362
18Fr x 2.3cm	Low profile balloon gastrostomy button kit	1	GBX1823	FAN85403
18Fr x 2.5cm	Low profile balloon gastrostomy button kit	1	GBX1825	FAN85404
18Fr x 2.7cm	Low profile balloon gastrostomy button kit	1	GBX1827	FAN85376
18Fr x 3.0cm	Low profile balloon gastrostomy button kit	1	GBX1830	FAN85363
18Fr x 3.5cm	Low profile balloon gastrostomy button kit	1	GBX1835	FAN85377
18Fr x 4.0cm	Low profile balloon gastrostomy button kit	1	GBX1840	FAN85345
18Fr x 4.5cm	Low profile balloon gastrostomy button kit	1	GBX1845	FAN85344
20Fr x 0.8cm	Low profile balloon gastrostomy button kit	1	GBX2008	FAN85405
20Fr x 1.0cm	Low profile balloon gastrostomy button kit	1	GBX2010	FAN85378


Nutricare® low profile balloon gastrostomy button kit codes

Size	Description	UOI	Code	NHSSC
20Fr x 1.2cm	Low profile balloon gastrostomy button kit	1	GBX2012	FAN85402
20Fr x 1.5cm	Low profile balloon gastrostomy button kit	1	GBX2015	FAN85373
20Fr x 1.7cm	Low profile balloon gastrostomy button kit	1	GBX2017	FAN85356
20Fr x 2.0cm	Low profile balloon gastrostomy button kit	1	GBX2020	FAN85401
20Fr x 2.3cm	Low profile balloon gastrostomy button kit	1	GBX2023	FAN85374
20Fr x 2.5cm	Low profile balloon gastrostomy button kit	1	GBX2025	FAN85365
20Fr x 2.7cm	Low profile balloon gastrostomy button kit	1	GBX2027	FAN85364
20Fr x 3.0cm	Low profile balloon gastrostomy button kit	1	GBX2030	FAN85381
20Fr x 3.5cm	Low profile balloon gastrostomy button kit	1	GBX2035	FAN85380
20Fr x 4.0cm	Low profile balloon gastrostomy button kit	1	GBX2040	FAN85382
20Fr x 4.5cm	Low profile balloon gastrostomy button kit	1	GBX2045	FAN85413
24Fr x 0.8cm	Low profile balloon gastrostomy button kit	1	GBX2408	FAN85395
24Fr x 1.0cm	Low profile balloon gastrostomy button kit	1	GBX2410	FAN85387
24Fr x 1.2cm	Low profile balloon gastrostomy button kit	1	GBX2412	FAN85368
24Fr x 1.5cm	Low profile balloon gastrostomy button kit	1	GBX2415	FAN85407
24Fr x 1.7cm	Low profile balloon gastrostomy button kit	1	GBX2417	FAN85386
24Fr x 2.0cm	Low profile balloon gastrostomy button kit	1	GBX2420	FAN85394
24Fr x 2.3cm	Low profile balloon gastrostomy button kit	1	GBX2423	FAN85371
24Fr x 2.5cm	Low profile balloon gastrostomy button kit	1	GBX2425	FAN85393
24Fr x 2.7cm	Low profile balloon gastrostomy button kit	1	GBX2427	FAN85412
24Fr x 3.0cm	Low profile balloon gastrostomy button kit	1	GBX2430	FAN85383
24Fr x 3.5cm	Low profile balloon gastrostomy button kit	1	GBX2435	FAN85347
24Fr x 4.0cm	Low profile balloon gastrostomy button kit	1	GBX2440	FAN85348
24Fr x 4.5cm	Low profile balloon gastrostomy button kit	1	GBX2445	FAN85406

Nutricare® feed sets for Nutricare low profile balloon gastrostomy button

The Nutricare feed set is an enteral feed set which attaches to the feeding/medicine port of the Nutricare low profile balloon gastrostomy button (Nutricare button). The Nutricare feed set is a non-sterile, single patient use device, made from silicone. All feed sets are fitted with an ISO 80369-3 compliant ENFit® connector.

Twice the usable life, compared with market leading alternatives*



Features and benefits

- **Silicone tubing**, not PVC, like many other brands of feed set
- **100% increase in use lifetime** compared to other feed sets on the market
- **Aids sustainability** targets and reduces plastic waste
- **Bolus and right angled configurations**



Size	Description	UOI	Code	NHSSC
280mm	Right angle feed set with ENFit connector and cap	5	GBF280R	FAE85180
280mm	Straight feed set with ENFit connector and cap	5	GBF280S	FAE85178
600mm	Right angle feed set with ENFit connector and cap	5	GBF600R	FAE85179
600mm	Straight feed set with ENFit connector and cap	5	GBF600S	FAE85181
90mm	Straight feed set with ENFit connector and cap	5	GBF90MS	FAE85182

*Based on internal testing. Actual duration may vary depending on clinical use conditions.

Nutricare® stoma measuring device

The Nutricare stoma measuring device is a sterile measuring probe made of high-quality medical silicone in size 12Fr. The Nutricare stoma measuring device determines the length of an already existing and healed stoma tract with an accuracy of ± 2 mm, in order to support the selection of the appropriate length of the gastrostomy button.



Features and benefits

- **12Fr size** used to measure the length of all stoma tract widths from 12Fr upwards
- **Medical grade** silicone tubing
- **Intuitive design**
- **Individually packed**
- **Free of latex**, PVC , DEHP and Bisphenol A

Size	Description	UOI	Code	NHSSC
12Fr	Nutricare stoma measuring device	4	NSMD	GPC85017

GBUK Home Hub

Delivery direct to your patient's door


Benefits of GBUK Home Hub:

- Guaranteed supply to your patients when they need it
- Control of ordering, re-ordering and amendments
- Quick and easy to set up or amend deliveries
- Scheduled timely deliveries according to patient needs
- Dedicated customer service support
- Local enteral specialist support

Patient benefits:


- Delivery guaranteed to their home address
- Call off and emergency orders available
- Never run out of product
- Customer helpline
- Support literature

Follow our simple steps:



GBUK Home Hub

Order online at www.gbukgroup.com using your smartphone, laptop or tablet.



GBUK Group Ltd.

Woodland House, Blackwood Hall Business Park,
North Duffield, Selby, North Yorkshire,
YO8 5DD, UK

T +44 (0) 1757 288 587

E info@gbukgroup.com

W www.gbukgroup.com

