Woodland House
Blackwood Hall Business Park
North Duffield, Selby, YO8 5DD

Compatibility statement

Glove Stand Aid Sling (GXP) and Comfort Value Stand Assist Sling (CVSA).

T

The GXP and CVSA are loop-attachment stand assist slings, designed for use with powered stand aids that are equipped to accept loop attachments. The slings attach to the stand aid via the designated loop attachment points.



* All slings, including their design features and intended application, are subject to an individual assessment to determine their suitability and safety for the person being transferred. This assessment includes consideration that the sling design provides adequate support, positioning, and transfer safety.

Registered office GBUK Group Ltd, Woodland House, Blackwood Hall Business Park, North Duffield, Selby, North Yorkshire, YO8 5DD. VAT registration number 228542895. A company registered in England & Wales with company number 09458917.