Woodland House
Blackwood Hall Business Park
North Duffield, Selby, YO8 5DD

Compatibility Statement

Comfort Value Transport sling.

The Comfort Value™ Transport Sling is designed for use with powered stand aids featuring a 4-point attachment system. The diagram (A) below provides an example of a compatible powered stand aid and highlights the attachment configuration required for the safe operation of the sling.

As a loop-attachment sling, the Comfort Value™ Transport Sling is compatible with powered stand aids that are designed to accept loop attachments.

Product Code- CVTR

Diagram A



* All slings, including their design features and intended application, are subject to an individual assessment to determine their suitability and safety for the person being transferred. This assessment includes consideration of sling appropriateness for individuals with amputations, including single above-knee and bilateral above-knee amputees, ensuring the sling design provides adequate support, positioning, and transfer safety.

Registered office GBUK Group Ltd, Woodland House, Blackwood Hall Business Park, North Duffield, Selby, North Yorkshire, YO8 5DD. VAT registration number 228542895. A company registered in England & Wales with company number 09458917.